



Azeer Sleep Diagnostic Devices



MiniScreen PRO

Description and instructions for use

HOSPITAL

HEMOCARE

DIAGNOSTICS

Table of Contents

1. Introduction.....	2
1.1. Assigned Purpose:.....	2
1.2. Contraindications.....	2
1.3. Accessories.....	2
2. Information.....	3
2.1. Safety-related information.....	3
2.2. General information.....	4
2.3. Technical inspection.....	4
2.4. Liability for functioning or damage.....	4
3. Recording measurement data using the MiniScreen Pro unit.....	5
3.1. Sensors.....	5
3.1.1. Sensors for flow and snoring noises.....	5
3.1.2. Thermistor sensor for breathing.....	5
3.1.3. Sensor for PAP (pressure).....	5
3.1.4. Sensor for oxygen saturation and pulse frequency.....	5
3.1.5. Sensor for thorax and abdominal movement.....	5
3.1.6. Sensor for body position.....	5
3.1.7. Sensor for leg movement.....	6
3.1.8. Sensor for ECG.....	6
3.1.9. Sensor for Snoring (external).....	6
3.1.10. Sensors for EEG / EOG / EMG.....	6
3.2. Preparing for measurement (Ambulatory).....	6
3.3. Attaching the MiniScreen unit and starting the measurement.....	7
3.4. Online Measurements (Inpatient).....	7
4. Service and Maintaining the device.....	7
4.1. Charging the battery.....	7
4.2. Cleaning instructions.....	8
4.3. Maintenance.....	9
4.4. Transport and storage conditions.....	9
5. Evaluation software for PC.....	9
5.1. Installing the software.....	9
5.2. Selecting a printer.....	9
5.3. Interfaces.....	9
5.4. Viewer Software "Sleep Xpert App".....	9
6. Troubleshooting.....	9
7. Ordering informations.....	11
8. Technical specifications.....	12
9. Used Symbols.....	13
10. EMC Information.....	14
11. Index.....	17

REV 2021-06-15

MiniScreen PRO

1. Introduction

1.1. Assigned Purpose:

The MiniScreen PRO Sleep Diagnosis System is a recording system for use in professional health care and home health care environment. It serves as a differentiated diagnosis for the sleep apnea syndrome and enables continuous recording of up to 40 channels for at least 15 hours without loss of data. It records signals for the following physiological variables:

- Flow
- Oxygen saturation SpO₂
- Pulse frequency
- Pulse wave
- Body position
- Snoring sounds (internal microphone)
- PAP (Positive Airway Pressure)
- Ambient Light
- Thorax effort
- Abdomen effort

optional:

- Thermistor
- 6 x ECG
- Central heart Frequency
- PTT (Pulse Transit Time)
- Systolic Blood Pressure Trend
- Snoring sound (external microphone)
- 2 x EMG (Leg movement)
- 3 x EMG (Chin)
- 6 x EEG
- 2 x EOG
- 8 x Analogue channel

1.2. Contraindications

There are no absolute or relative contraindications to using the device. In the following cases, the device must be used under the supervision of qualified medical personnel:

- in patients with acute life-threatening illnesses.
- in patients with acute, serious infections.
- in mentally confused patients.
- in infants and children

1.3. Accessories

The MiniScreen Pro system as a whole (ME system) contains three components:

- The microprocessor controlled MiniScreen Pro recording device
- PC software for measurement data presentation and PC analysis
- Optional: Analogue box for 8 analogue inputs incl. connection cable

It is also equipped with the following applied parts:

- SpO₂ finger sensor with cable for obtaining pulse frequency and oxygen saturation values

- Velcro strap for securing the finger sensor to the wrist
- Flow prongs for obtaining the respiratory signal. The 20 cm long adapter hose is identified by a blue ring
- Adapter hose for measurement during PAP Therapy
- Flexible carry belt with integrated pressure transducers (thorax effort) to fasten the MiniScreen Pro unit to the patient
- Flexible belt with integrated pressure transducer to record abdominal effort
- Shoulder bag for storing the MiniScreen Pro unit

Other accessories:

- USB Interface cable for data transfer between the MiniScreen PRO and PC
- Battery recharger
- Transport case

optional applied parts:

- ECG electrode for recording ECG signals
- Thermistor for additional breathing detection
- External snoring microphone

Other optional accessories:

- Splitter box for connecting the Leg sensors
- Splitter box for connecting the EEG/EOG/EMG signals
- USB isolator for online measurement
- Viewer "Sleep Xpert App" for tablet.

The PC, the charger and the Analogue box are not suitable for a domestic or patient environment due to lack of drip protection!

The data can be displayed and analysed on a standard commercial PC. The measurement curves and the analysis results can be printed on all commonly available printers, including dot-matrix, laser or ink jet printers.

The MSV (MiniScreen Viewer) analysis program has the following minimum configuration requirements

- PC with OS Windows 7 Professional, Windows Server 2012 R2 or higher
- 8 GB of RAM (32 GB recommended)
- 1 GB of disk space on the hard drive
- Mouse
- Dedicated graphics card with a resolution of at least 1920 x 1200 and 32 bit colour resolution
- Free USB port
- Printer with Windows driver

2. Information

2.1. Safety-related information

Observe the instruction manual:

Every use of the device requires exact knowledge and observance of this instruction manual. The device is only intended for the use described herein.

No alarm function available!

The device is not suitable for continuous monitoring of vital or physiological functions (for example, intensive monitoring, monitoring operation) since there is no SpO2 alarm. No direct data analysis takes place in the device.

No SIDS Monitoring:

The device is not suitable for use as a SIDS monitor (SIDS: Sudden Infants Death Syndrome, sudden child death).

Patient instruction:

Patient instruction must be provided by a healthcare professional trained in the device. An attached short manual does not replace the manual instructions or the warning of possible hazards.

Separation from the supply network:

To isolate the device from the supply network, the plug of the power supply must be disconnected.

Do not open the device!

Warning:

Additional devices that are connected to medical electrical (ME) devices must be able to meet their IEC or ISO standards. Anyone connecting additional devices to ME devices is a system configurator and is responsible for the ME system being in compliance with the normative requirements (for example, IEC 60601-1).

Warning:

For online measuring you must use an electrical or optical isolator for connection of the device to the PC (available from your home care provider). A connection of the device to the PC without electrical or optical isolator is only allowed if all patient connections have been removed first. Only Physicians or trained staff is allowed to perform online measurements.

Warning:

Magnetic and electrical fields can impact on the functioning of the device. Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally. When operating the

device, ensure that all third-party devices being operated in the vicinity are meeting their relevant EMC requirements. X-Ray equipment, HF surgical devices, tomography devices etc. can interfere with other devices as they may emit higher levels of electromagnetic interference.

Warning:

The device has no defibrillator protected application parts! Disconnect the device before defibrillation! A direct application of the device to the heart (esp. ECG!) is not allowed!

While the device is being worn, no invasive or intracorporal measurements and interventions (e.g., electro- or HF-surgical devices) may take place.

Warning:

Using multiple devices for patient screening could result in excessive leakage current above the allowed value!

Warning:

When using the device on patients with active implants, such as Cardiac pacemakers, peripheral nerve stimulation devices, tongue pacemaker etc., medical professionals should be alert to possible interference with the device or the implant.

Warning:

Applying the device to more than one patient at any one time is not allowed.

Warning:

The device must not be operated in vehicles and aircraft.

Warning:

When applying sensors to the patient, please ensure that no electrical line of the device is in contact with any other electrical conductive parts including ground.

Warning:

Avoid placing this instrument in direct sunlight or in close proximity to intense heat. Prevent also the contact with dust, lint, dirt, moisture and liquids.

Warning:

Children or incompetent persons should not use this device unattended without first having obtained instructions for the safe use of the device. Ensure that infants, children and animals cannot touch the device unattended

Warning:

Caution must be taken to ensure that cables do not encircle the patients neck. Special attention is needed in the case of children.

 Warning:

Before the battery is charged, the power adapter must be checked for external damage!

Never use the power adapter with damaged housing or cable!

 Warning:

Before every application, the housing of the device as well as the cables and sensors must be checked for external damage!

Never use the device with damaged housing, cables or sensors!

 Warning:

The device and PC software are not intended for automatic diagnosis only. The measurement data must be always reviewed and evaluated manually by a qualified physician.

2.2. General information

This instruction manual should be regarded as a component of the device. It should be kept on hand somewhere near the device at all times. Reading and understanding the instruction manual is a prerequisite for proper use and correct handling of the device in order to maintain the safety of the patient and operator.

The guarantee is valid for a period of 24 months for the device and 6 months for the equipment from the date of sale.

Only accessories which are listed in this instruction manual should be used with the device. We cannot guarantee the safe operation and function of the device if unknown proprietary accessories / consumables are used that have not been tested. (E.g. Patient leads, Sensors, Consumables, Memory cards etc.)

Damages resulting from the use of third-party accessories or consumables shall render this warranty void.

The manufacturer will only assume responsibility for the device in terms of safety, reliability and function if:

- 1. assembly, add-ons, reinstallation, changes and repairs are carried out by the manufacturer or a qualified agent authorized by the manufacturer to do so;**
- 2. The device is used in accordance with the instruction manual.**

All printed material relates to the model of the device and the safety regulations at the time of printing. All devices, switches, processors, software programs and names contained herein are subject to copyright law.

LM shall only be liable for the malfunction of the device and its software if used in the normal operating conditions in accordance to this manual.

If a PC is included in delivery third party software is not allowed to be installed on this PC.

Medical devices must only be operated by a qualified experienced person(s) to ensure correct handling of such devices.

The operator must read and understand the user manual to operate the device correctly.

The operator must check the functionality of the device before each use to ensure it is in sound condition and in good working order.

Functional testing of the device must be carried out at regular intervals. It is recommended that this is conducted once a month.

Upon reaching the end of its service life, the device and its accessories should be disposed of in accordance with the WEEE Directives or relevant Electronic Disposal protocol.

For further details please contact LM.

2.3. Technical inspection

Only devices that are regularly checked and maintained are deemed safe to use. It is recommended that units are subject to a test protocol and the battery is replaced every 24 months.

For more details on servicing and service contracts, please contact the LM Customer Service Department or your service provider. A functional tester can not be used to assess the accuracy of the finger sensor or pulse oximeter of the device.

The following tests need to be carried out by the operator before each measurement:

- Visual inspection of the device and accessories for obvious damage which could result in mechanical failure.
- Testing all hose connections are air tight.
- Testing the thorax and abdominal sensors for leaks.
- Checking the display LEDs.
- Testing the pulse and SpO2 finger sensors and carrying out a plausibility test (timed pulse test using).
- Testing the battery capacity.

2.4. Liability for functioning or damage

In the event of improper use or repair by the owner or operator, the liability for the device will be transferred to those parties. LM shall not be liable for damages arising from non-observance of this instruction manual. The actual guarantee and warranty conditions in the LM terms of sale/delivery are not extendable.

3. Recording measurement data using the MiniScreen Pro unit

3.1. Sensors

3.1.1. Sensors for flow and snoring noises

The detection of patient flow is determined by the use of a nasal flow prong that transmits the pressure signal to an internal pressure transducer of the device. The flow prong can be easily applied by the user and will have no adverse affect on the quality of sleep. The use of a flow prong also eradicates the problem where adhesive sensors cannot be used. (i.e.: patients with beards). In addition the flow prong is an economical choice of sensor as it is commonly available.

No additional sensor is required for respiratory and snoring detection. Respiratory sounds are transmitted to the MiniScreen Pro via the flow prong. The built in microphone detects the sounds which are then electronically analysed. Therefore additional microphones are not required.

Due to the high sensitivity of the internal pressure transducer, the MiniScreen Pro is capable of measuring and recording ultra fine pressure differences associated with patients who mouth breathe.

In order to connect the flow prong, ensure the connector of the prong attaches into the **blue** connector of the device.

The flow prong should be used in accordance with the manufacturer's instructions.

Warning:

The flow prong is intended for single patient use. No attempt should be made to clean or use the flow prong the more than once. Multiple use can result in cross infection!

3.1.2. Thermistor sensor for breathing

Additional to the flow prong a Thermistor can be used to measure the breathing. Please also observe the instructions accompanying the sensor.

The **brown** connector of the Thermistor needs to be attached to the **brown** connector on the MiniScreen.

3.1.3. Sensor for PAP (pressure)

When performing measurements on a patient undergoing PAP therapy, an adapter hose is used instead of the flow prong.

The connector of the PAP adapter needs to be attached to the **blue** connector on the MiniScreen.

3.1.4. Sensor for oxygen saturation and pulse frequency

A pulse oximeter has been integrated into the MiniScreen Pro for the purpose of measuring oxygen saturation and pulse frequency. A failure of the finger sensor or a missing data update by the pulse

oximeter is indicated to the user by the red LED on the device. At the same time the values for SpO2 and pulse are fed to 0.

Please ensure that blood circulation at the point of measurement is not affected by the manner in which the sensor has been secured. Excessive pressure should not be exerted on the finger, especially if the temperature exceeds 41 ° C.

When using a finger sensor, please use the Velcro wrist/armband to secure the sensor and cable. All nail polish (coloured or clear) must be removed from the measuring finger to obtain usable data.

To minimize interference (e.g., motion artifacts), the pulse oximetry values are filtered by digital data processing. As a result of data averaging and signal processing, slight delays in the display of the pulse oximetric values occur. The internal pulse oximeter works with a time horizon of 4 seconds. By additionally considering the change tendency, the minimum value of saturation at the end of an apnea is correctly reproduced. The data is updated at every beat, so there are no measurable delays due to data refresh and transmission. Please read and understand the instructions accompanying the finger sensor.

3.1.5. Sensor for thorax and abdominal movement

The sensors used for recording thorax and abdominal movement consists of small rubber pressure pads that are connected to the MiniScreen Pro via measuring hoses. The sensor for recording thorax movement comprises of two pressure pads whereas the sensor for abdominal movement contains only one. The sensors are inserted into the pockets of the elasticated effort belts. The thorax belt is applied at the height of the sternum, the abdominal belt in the stomach region.

For hygiene reasons and to avoid allergic reactions the belts should be **worn over clothes**.

The **red** connector of the thorax sensor needs to be attached to the **red** connector and the **black** connector of the abdominal sensor needs to be attached to the **black** connector of the device.

The effort belts are adjustable in size due to the Velcro fasteners, and should be suitable for most patients. These belts are however also available in special sizes.

Note:

Overstretching of the belts may result in poor signals and possible loss of data.

3.1.6. Sensor for body position

The position sensor integrated in the MiniScreen Pro displays information related to the patient's body position during the study.

For accurate determination of body position, ensure the MiniScreen Pro is applied to the patient correctly.

MiniScreen PRO

The body positions detected are: Supine, Prone, Left side, Right side, Upright.

3.1.7. Sensor for leg movement

To diagnose restless or periodic leg movements (restless leg), the MiniScreen Pro can be equipped for each leg with a sensor and internal recording software (optional). The device then enables continuous recording of leg movements and associated analysis results in the report.

Attaching the leg sensor:

All commonly available types of electrodes with 1.5mm safety plugs can be used for both bipolar EMG channels. To connect the electrodes to the device the splitter box is used. The **yellow** marked connecting plug of the splitter box needs to be plugged into the **yellow** socket on the MiniScreen.


For application and cleaning of the electrodes please refer to the manual of the manufacturer.

For long time electrodes a very careful preparation of the skin is necessary!

3.1.8. Sensor for ECG

The supplied ECG cable (optional) is compatible with all commonly available adhesive electrodes and can be used as consumables.

Attaching the ECG electrodes:



The **red** marked connecting plug of the EEG electrode cable needs to be plugged into the **red** socket on the MiniScreen Pro.

Important: The ECG is not suitable for a differential cardiac diagnosis!

3.1.9. Sensor for Snoring (external)

Additionally to the internal snoring microphone an external microphone can be used to obtain the snoring sound. For adapting the snoring microphone please refer to the manual of the sensor.

The **blue** marked connecting plug of the microphone needs to be plugged into the **blue** socket on the MiniScreen.

3.1.10. Sensors for EEG / EOG / EMG

For the neurological channels (EEG, EOG, EMG) all commonly available types of electrodes with 1.5mm safety plugs can be used. To connect the electrodes to the device the splitter box is required. The **green** plug of the splitter box has to be connected to the **green** socket of the device.

For application and cleaning of the electrodes please refer to the manufacturers instructions.

Note: In order to guarantee good adhesion of the electrodes and strong frequency signals, clean the desired skin area with a sterile solution and apply a good quality conductive cream

3.2. Preparing for measurement (Ambulatory)

In preparation of an ambulatory study where the patient will use the device at home, please ensure the following:

1. Charge the battery: see (Page 7).
2. Connect the MiniScreen Pro with the USB interface cable.
3. Start the MiniScreen software on the PC.
4. Fill out the fields relating to personal data of the patient in the menu item "Record / Initialize device (Offline)" and start the transfer.
5. Switch on the MiniScreen Pro by sliding the switch to the right. For testing purposes, both LEDs will light up initially. The red LED will then go off. The green LED on the MiniScreen Pro will remain lit for the duration of the initialization process.
6. In the next step of the Initialization process, you have the option to choose Timer Controlled or Manual recording.

Timer controlled recording start:

A dialogue window will be displayed in order for the operator to determine the start time for recording. The MiniScreen Pro needs to remain **switched on** after this time!

The unit will now rest in stand-by mode until the selected Start Time of the recording. The unit will then automatically switch on and start recording data. **Never switch off the MiniScreen Pro sliding switch during use!**

Manual recording start:

Once the unit has been successfully initialized both LEDs will flash confirming ready for use. The MiniScreen Pro unit now needs to be switched off sliding the switch to the left. After fitting the device the patient will need to switch the unit on before going to bed by sliding the switch to the right.

7. The MiniScreen Pro is now ready for use. The unit can now be disconnected from the PC.
Note: remove cable by pulling on the plug - not on the cable.

3.3. Attaching the MiniScreen unit and starting the measurement

When apply the sensors, make sure that they are applied in such a way that strangulation is not possible due to their length (for example flow prongs, finger sensor). In order to ensure reliable recording of data, you need to observe the following points when attaching the MiniScreen Pro and sensors. Before ambulatory measurements these steps should be demonstrated and practised with the patient beforehand in the clinic:

1. If applicable, attach the ECG electrodes.
2. Secure the MiniScreen Pro to the thorax using the elasticated thorax belt (with the **two** built in pockets). Apply the belt whilst in the standing position and exhale. A correctly fitted belt should not slip up or down once fastened. To record abdominal breathing, attach the abdominal belt (with the **single** built-in pocket) around the stomach in the same way. Always ensure the belts are fitted over the night clothes and not directly onto bare skin.
3. Check the position of the pressure transducers (black rubber pads). The sensors should be fully inserted into the built-in pockets sewn onto the belt. **Never pull on the tubing connected to the pressure pads.**
4. Apply the flow prong and/or respiratory thermistor to the nose or fix the adapter hose to the oxygen mask.
5. If applicable, attach the leg sensor, the microphone and the neurological splitter box.
6. Apply the finger sensor (SpO2) to a finger and secure the cable to the wrist using the Velcro wrist band. Whilst doing so, ensure that you do not apply pressure to the finger thereby disturbing the blood flow.
7. If no automatic start was predefined the MiniScreen Pro has to be switched on using the slider (towards the right). Otherwise the MiniScreen will automatically switch itself on at the predetermined time. Both LEDs light up briefly for testing purposes.
If all sensors have been connected correctly, the red LED will go out. The green LED blinks at regular 4 second intervals.
The red LED signals a faulty pulse signal. Check the finger sensor and cable.

Once the patient has been educated on fitting the device, the device can then be packed into its hard carry case for the patient to take home. The patient can refer to the quick reference guide when fitting the device at home.

When going to bed the patient needs to:

- Attach the MiniScreen Pro using the elasticated thorax belt.

- Attach the flow prong and finger sensor.
- If applicable, attach the abdominal belt, leg sensor, ECG and EEG electrodes.
- Check that the sensors and their connectors are positioned correctly.
- If manual recording start was selected: Switch on the MiniScreen Pro by sliding the switch to the right.

The next morning:

- Stop the measurement by switching off the MiniScreen Pro.
- Detach the MiniScreen Pro and sensors and return all parts to the case. Cleaning of the sensors will be carried out by the technical staff at the clinic or hospital.
- Bring the case back to the clinic or hospital.
- For automatic evaluation the measurement can be transferred to the PC using the device software. Select the "Record" tab followed by "Read measurement from device".

3.4. Online Measurements (Inpatient)

Please adhere to the following steps when performing an online measurement:

1. Ensure the battery is fully charged (see Page 7).
2. Apply the device and connect all sensors to the patient (see Page 7)
3. Connect the MiniScreen unit (switched off) to your PC using the USB interface cable.
4. Start the MiniScreen software on your PC.
5. Complete the Patient data fields in the menu item "Record / Start recording (online)" and start the transfer.
6. Switch on the MiniScreen Pro unit by sliding the switch to the right position. The device can be switched on also direct after attaching the device to the patient.
7. Perform impedance test and bio signal calibration (see the online help of the software).
8. Once the test is complete, (the following morning) click the stop button to end the measurement.
9. Switch the MiniScreen Pro device off by sliding the switch to the left position.

4. Service and Maintaining the device

4.1. Charging the battery

The MiniScreen Pro unit is equipped with a special fast-charge Li-Ion battery. The battery charging adaptor included with your device is specially designed for this type of battery and should only be used with your MiniScreen Pro. It is recommended to charge the battery completely after each measurement.

Caution: It is not allowed to use a different battery charger than the one provided!

MiniScreen PRO

Battery charging procedure:

- Connect the battery charger adaptor to the MiniScreen Pro
- Plug the charger into a power outlet
- The charger LED will illuminate white. As soon as the battery is charged, the LED will illuminate blue. The device can remain connected to the charging adaptor for a long period of time without suffering any damage.

4.2. Cleaning instructions

General

As with any medical device, when using the MiniScreen PRO polysomnographic system, certain hygienic procedures are required for safe reuse on the patient. Reusable products must have safe disinfectability to exclude any risk of infection from subsequent users / patients. The regulations of the Medical Devices Act stipulate that for such medical devices disinfection activities with procedures according to the RKI guideline are to be carried out in analogy to the methods of surface or instrument disinfection. It is recommended to create a hygiene plan for the MiniScreen PRO product. Sterilization of the products listed below is not necessary.

The MiniScreen PRO polysomnographic system is a medical device that is in the patient's environment when used as intended, with direct contact with the hands of the staff and the patient. These products are therefore classified as so-called "non-critical medical products" in accordance with the guidelines for hospital hygiene and infection prevention.

The recommendations are based on the following provisions:

Lists of tested and approved disinfectants and procedures by VAH, IHO, DGHM, RKI

As a disinfectant of the applied parts a product should be selected for rapid disinfection with short exposure times on an alcohol-free basis, which is listed in VAH, IHO or DGHM (for example MikroZid® sensitive). The information given by the manufacturer regarding dosage and time of use should be considered.

Periods

Applied part	Cycle		Cleaning		Disinfection	
	After each use	After patient change	By hand	Washing machine	Spray disinfection	Wipe disinfection
Basic device MiniScreen PRO		•	•			•
Carrying case, outside		•	•			•
Carrying case, inside		•	•		•	
Shoulder bag cloth bag with		•	•			•

shoulder strap						
Connection cable MiniScreen PRO	•		•			•
SpO2 finger sensor	•		•			•
Velcro wrist strap for finger sensor	•		•			•
Thorax- /Abdomen sensor		•	•			•
Thorax- /Abdomen strap		•		•		•
Electrode extension cable	•		•			•
Thermistor	•		•			•
LEMG Adapter cable	•		•			•
Snoring microphone	•		•			•
ECG electrode cable	•		•			•
Splitter box for EEG/EOG/EMG	•		•			•
Flow prong	Disposables (Single use)					
CPAP Adapter	Disposables (Single use)					
Electrodes with / without integrated connection cable (except gold cup electrodes)	Disposables (Single use)					

Cleaning

It is recommended that all components of the MiniScreen PRO be cleaned prior to disinfection to remove coarse dirt, adhesive residue, odours and the like. All components of the MiniScreen PRO can be wiped with a slightly damp cloth and a light detergent. Allow all components to air dry completely before reinstalling.

⚠ Warning:

No liquid should enter the MiniScreen PRO or plug connections. The sensors must not be immersed in liquid. Before cleaning, close the pressure connections of the device with caps!

Disconnect the charging and data cable from the device before each cleaning.

In the washing machine

The thorax and abdominal straps can be washed in the washing machine at 60° C (use a laundry net!). Let the straps dry in the air.

Attention:

- Do not use a dryer
- Remove thoracic and abdominal sensor beforehand

Disinfection

It is recommended to subject all components of the MiniScreen PRO to a disinfection procedure. The periods can be found in the table. In general, a product should be used for rapid disinfection with short exposure times on an alcohol-free basis. Observe the instructions for use of the disinfectant used.

- Wipe and spray disinfection with MikroZid® sensitive liquid
- Dive disinfection with gigasept® FF (new)

Disposables

Disposables must not be used multiple times or on different patients. Follow the guidelines for waste disposal in hospitals.

4.3. Maintenance

We recommend the device including all accessories is serviced every 2 years. The servicing should only be carried out by the manufacturer or by an authorised agent. Calibration of the pressure channel should be carried out once at the first time when the device is used. In order to maintain basic safety and essential performance, no special measures need to be taken with regard to EMC.

4.4. Transport and storage conditions

The ambient temperatures for transport and storage lie between:

-25 ° C and + 5 ° C,
+ 5 ° C to + 35 ° C with a relative humidity of up to 90%, without condensation,
> 35 ° C to 70 ° C at a water vapor pressure of up to 50 hPa.

Condensation must always be avoided.

5. Evaluation software for PC

5.1. Installing the software

To install the MiniScreen Pro software, place the CD supplied with your device into the CD-ROM drive of your computer. The program will automatically run. If the program does not run automatically open the Explorer option on your operating system and select the CD-ROM drive. Once selected double click the Setup.exe file from the software.

You can check your computer and network (if applicable) for MiniScreen software versions that may have been previously installed. To ensure the correct up to date version is selected, open the software installed and click the Search button. The program will then automatically find and use the latest version installed.

After selecting the target drive the installation program will now automatically create the necessary directories and copy all the required files to the hard drive. The installation will automatically insert a new group containing an entry for the new software in the Programs folder of the Start Menu. In addition a shortcut on your Windows desktop will be created.

After successful installation, remove the CD and store it in a safe place.

5.2. Selecting a printer

The measurement curves and reports can be printed on any commonly available type of printer. The correct printer driver can be selected using the Windows Control Panel application.

5.3. Interfaces

The MiniScreen PRO software has analogue and digital interfaces for import of external measurement data. The data can be imported by communication with a device (e.g. network, analogue input) or by data file. The imported data can be displayed in the

raw data screen and in the report together with the data of the MiniScreen or as stand alone data.

5.4. Viewer Software "Sleep Xpert App"

During an online measurement, the measurement curves can be displayed on a tablet (Windows or IOS) optionally with the help of the software "Sleep Xpert App". The tablet only serves to visualize the measured curves. There are no data stored and no findings or evaluation can be made.

6. Troubleshooting

Channels (e.g. pressure) are missing in the display.

The channels are deactivated and therefore do not appear during "Test" and "Record".

Check which channels are active for recording using the menu item Options / Edit Channel sets.

Channels are missing after loading a measurement.

They were not recorded or have been hidden for measurement data display.

Flow signal is missing during recording or follows the edge of the range.

The flow signal follows a straight flat line in the middle of the channel.

Check the flow prong on the patient and the connection on the MiniScreen unit. When measuring during PAP therapy, check the connection of the PAP adapter hose on the oxygen mask and on the MiniScreen Pro unit.

Signal amplitudes are very small or do not appear.

Check the corresponding sensors on the patient and their connections to the MiniScreen Pro unit. Whilst doing so, check that the tube hoses and the black pressure pads of the thorax sensors are intact. A leaking sensor can cause inaccurate readings in the thorax channel. The pressure pads should be attached to the patient securely without over tightening the belt and the tube hoses should be routed in such a way that they cannot kink.

Check all cables and connectors on the MiniScreen Pro unit and the PC.

Pulse oximeters not responding.

The channels for oxygen saturation and pulse frequency are registering 50% and 30 P/min respectively despite the patient being connected.

First check that the finger sensor is seated correctly on the patient and remove any nail polish that may be present. A small red light should be present in the sensor when the finger is inserted. If the light doesn't switch on when the sensor is applied to the finger, then you need to check the connection at the MiniScreen Pro unit and any intermediate extensions.

Printout not working.

The printer prints characters on the page but in no apparent format.

The wrong printer or printer driver has been installed.

MiniScreen PRO

The printer is not responding to the print command.

Check the printer cable and the connections to the printer and PC. The printer should be switched on and operational, i.e. the control LED's on the printer should be lit.

Cannot establish connection to the MiniScreen.

The USB cable is not connected correctly.

Check the connection of the USB cable to the MiniScreen Pro and the PC.

USB interface was deactivated

The USB interface can be activated in the MiniScreen Pro software. Select the "Options" tab followed by menu "Device Settings".

Battery is discharged.

Charge the battery correctly.

For clarification and in case of problems during the installation, maintenance or use, please refer to your local distribution partner, or directly to the manufacturer.

7. Ordering informations

Art.-No.	Spare parts / consumer material / Accessories
200-0312/01	Nasal Cannula Adult, 20 cm, open connection
200-0312/10	Nasal Cannula Adult, 50 cm, Luer-Lock
200-0155/22	CPAP adapter with 150cm tube
200-0313/01	Adapter tube, luerlock-male 20 cm for CPAP-Adapter 200-0155/22 + Nasal Cannula Adults, 20cm 200-0312/01
100-0031	ECG electrode fleece, 50 cm cable (Pack à 12)
100-0031/02	ECG electrode, children fleece (Pack à 12)
100-0031/01	ECG-electrodes for PSG wet gel, foam, push button (Pack à 25)
512 285	EEG-adhesive electrode PSG (Pack á 25)
100-0032/12	EEG-cup electrodes for PSG 100cm, different colours (Pack á 10)
100-0032/11	EEG-electrodes for PSG 100cm, white connector (Pack á 10)
100-0032/01	EEG-electrode for PSG - white (Pack á 10)
100-0031/10	Adhesive pads for electrode application (Pack á 50)
100-0051/01	Abrasive cleaning paste, SkinPrep 130g
100-0051	Abrasive cleaning paste, SkinPure 135 g
510001	Electrode Cream SAC2 tube 100g
150-0005	EEG electrode cream Ten20 114 gr., tube
150-0005/01	EEG electrode cream Ten20 114 gr., can
200-0316	Thorax sensor for adults MiniScreen Plus / Pro
200-0917	Strap thorax size S, 98 cm blue, MiniScreen Pro
200-0918	Strap thorax size M, 122cm black, MS Pro
200-0919	Strap thorax size L, 142 cm red, MS Pro
200-0320	Abdominal sensor for adults MS Plus / Pro
200-0921	Strap abdomen size S, 98 cm blue, MS Pro
200-0922	Strap abdomen size M, 122cm black, MS Pro
200-0923	Strap abdomen size L, 142 cm red, MS Pro
200-0351	Repair set for MiniScreen abdomen/thorax sensors
200-0311	Arm wrap for SpO2 sensor
200-0310	SpO2 softclip fingersensor latex free, MiniScreen
200-0952	ECG electrode cable red / yellow / black
200-0958	LEMG adapter for MiniScreen Pro
100-0119/03	extension cable LEMG 100 cm safety connector/clutch
200-0957	EEG / EOG / EMG adapter, MiniScreen Pro
200-0861	Flowsensor oral / nasal, white 1m, cannula mount
200-0930	Snore microphone, white, MiniScreen Pro
200-0910	Battery charger with lemo plug, MS Plus / Pro
200-0936	USB connecting cable, MiniScreen Plus / Pro
200-0912	Patientbag, MiniScreen Pro
200-0913	Transport case with inlet, MiniScreen Pro
160-0308/01	Adhesive patch for snore microphone (Pack á 50)

Art.-No.	Spare parts / consumer material / Accessories
160-0300/01	Skin marker set
160-0301/01	Löwenstein Medical Measuring tape
160-0302	Löwenstein Medical Scissor
160-0303	LM tongue depressor non-sterile (Pack á 100)
160-0304	LM cotton swab non-sterile (Pack á 100)
160-0306	ES compresses non-sterile 10x10cm (Pk á 100)
160-0307	Medipore fixation fleece hypoallergenic, 10 cm x 10 m
160-0308	Durapore plaster hypoallergenic, 2,5cm x 9,1m
160-0309/01	Disinfection towels can á 200 towels

Art.-No.	Spare Parts / Accessories for use with children
200-0312/31	Cannula for dynamic pressure measurement, children, 50cm
200-0368/05	Dynamic pressure nasal cannula infant, 30 cm
200-0368/06	Dynamic pressure nasal cannula neonate, 30 cm
200-0368/08	Dynamic pressure nasal cannula pediatric, 30 cm
200-0155/07	Adapter tube, 100 cm
200-0364	Sensor thorax children MiniScreen Plus/Pro
200-0965/03	Strap thorax size XS, 79 cm conn. for position sensor
200-0965/04	Strap thorax size XXS, 65 cm connector for position sensor
200-0965/05	Strap thorax size XXXS, 40 cm connector for position sensor
200-0366	Sensor abdomen children, MiniScreen Plus/Pro
200-0967	Strap abdomen, size XS, 79 cm, MiniScreen Pro
200-0967/01	Strap abdomen size XXS, 65 cm, MiniScreen Pro
200-0967/02	Strap abdomen size XXXS, 40cm, MS Pro
300-0201/01	SoftClip finger sensor, 45cm paediatric, 15-50 kg *in combination with 200-0370
300-0201/02	SoftFlex-sensor neo. 90 cm 1 - 4 kg, foot and finger *in combination with 200-0370
200-0310/08	Glue Sensor SpO2 neonatal, 1-3 kg, 50cm *in combination with 200-0370
200-0310/09	Glue Sensor SpO2 Infant, 3-15 kg, 50cm *in combination with 200-0370
200-0310/10	Glue Sensor SpO2 paediatric, 15-45 kg, 50cm *in combination with 200-0370
200-0370	Connection cable SPO2

MiniScreen PRO

8. Technical specifications

- Dimensions : 35mm x 75mm x 168mm (H x W x L, without bag)
- Weight : 260 g including storage battery, without bag
- Housing : metallized plastic (ABS, UL 94HB)
- Temperature range : +5°C...+40°C
- Moisture : 10% - 90%
- Atmospheric pressure: 70kPa - 106kPa
- Storage media : Internal flash memory
- Storage capacity : min. 24 hours
- Registered parameters :
 - Respiratory activity : Differential pressure measurement via flow prong (with adaptor also during PAP therapy), alternative or additional measurement by means of thermistor possible (option)
 - Thoracic effort : Differential pressure measurement by means of rubber cuffs built into chest strap
 - Abdominal effort : Differential pressure measurement by means of rubber cuffs built into abdomen strap
 - Breathing sounds : Phonometric transducer via flow prong
Alternative or additional measurement with external microphone possible (option)
 - SpO₂/Pulse : Built-in pulse oximeter, calibrated for functional oxygen saturation
SpO₂ measurement range: 80%-100% ± 2% SpO₂ 60%-79% ± 4% SpO₂
Pulse measurement range: 50 1/min - 150 1/min ± 2%
(Reference: electrical pulse simulator)
Finger sensor: special rubber-coated thimble finger sensor
 - Pulse wave : Plethysmogramm display; measurement via finger sensor
 - Position : Acceleration sensor for position recording (5 positions)
 - Ambient Light : Photometric measurement and light-Intensity display
 - PAP : Differential pressure measurement directly on oxygen mask
Measurement range: 0cmH₂O - 45cmH₂O ± 5%
 - Leg movement : Two leg sensors (EMG) for measurements of muscular actions, separate recording for left and right leg; connection via 1.5mm safety plug (option)
 - ECG : 6 channel lead via adhesive electrodes (option), suitable for patients <10kg
 - Central heart frequ.: Measurement range: 30 1/min - 200 1/min ± 2% (option)
 - PTT : Measurement range: 100 ms – 355 ms ± 4% (option)
 - Syst. blood pressure : Trend (option)
 - EEG : 6 channel lead via adhesive electrodes (option); connection via 1.5mm safety plug;
Impedance: 10 MΩ, frequency: 0,2 Hz - 45 Hz
 - EOG : 2 channel lead via adhesive electrodes for left and right eye (option); connection via 1.5mm safety plug
 - EMG : Bipolar measurements (3 x EMG) of muscular actions of the chin via adhesive electrodes (option)
 - External channels : External box with voltage input (RJ11; 0..2.5 V) for up to 8 external channels with galvanic separation and RJ11-jack
- Fault indicator : 2 LEDs on front of the device
- Power supply : Rechargeable Li-Ion storage battery 3.6 V with built-in Semiconductor safety
- Charger : Plug-in power supply with medical approval
- Output : USB interface with cable for data transmission
- Power consumption : Approx. 160mA
- Online operations : In online operation with a patient, an optical waveguide to the PC is indispensable (option)

MiniScreen PRO

9. Used Symbols

<i>Symbol</i>	<i>Meaning</i>	<i>Symbol</i>	<i>Meaning</i>
	Pay attention to the instructions for use!		Electrical and electronic devices may not be disposed of with domestic waste. Consumers are obliged by law to return electrical and electronic devices at the end of their service lives to the public collection points set up for this purpose or point of sale. Details to this are defined by the national law of the respective country. This symbol on the product, the instruction manual or the package indicates that a product is subject to these regulations. By recycling, reusing the materials or other forms of utilising old devices, you are making an important contribution to protecting our environment.
	Pay attention to the instructions for use!		
	Type BF		
	Protection class II		
	Best before date		
	Manufacturer	IP22	Device is protected against solid foreign bodies with diameter ≥ 12.5 mm. Device is protected against access with a finger. The device is protected against falling dripping water when the housing is tilted up to 15° .
	Manufacturing date		Once product! Neither suitable for a preparation nor for a multiple application.
	Medical device		

MiniScreen PRO

10. EMC Information

Medical electrical devices subjected to special preventive measures regarding to electromagnetic compatibility (EMC). They must be put in operation and installed according to the EMC information in this manual. Portable and mobile RF communications equipment (e.g. mobile phones) can influence medical electrical devices. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the MiniScreen PRO, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. For an intended use of the MiniScreen PRO, the specified equipment in this manual must be used. The application of other equipment can lead to a higher emission and a reduction of the immunity.

Guidelines and Manufacturer's Declaration – Electromagnetic Emissions		
The MiniScreen PRO is intended for use in an electromagnetic environment as specified below. The customer or user of the MiniScreen PRO should ensure that it is used in such an environment.		
Emission Test	Compliance	Electromagnetic Environment - Guidelines
RF – Emission according to CISPR 11	Group 1	The MiniScreen PRO uses RF energy only for its internal function. So its RF emission is very low and it's implausible that it cause any interferences to surrounding electronical devices.
RF – Emission according to CISPR 11	Class B	
Harmonic Emissions according to IEC 61000-3-2	Not applicable	The MiniScreen PRO is suitable for the use in all establishments, including domestic establishments and those, which are connected to the public power supply, which powers buildings for residential purposes.
Voltage Fluctuations / Flicker Emissions according to IEC 61000-3-3	Not applicable	

Table 1: Table 201 EN 60601-1-2, Electromagnetic Emissions

Guidelines and Manufacturer's Declaration – Electromagnetic Immunity		
The MiniScreen PRO is intended for the electromagnetic environment as specified below. The customer and user of the MiniScreen PRO should ensure that it is used in such an environment.		
Immunity Test	Test Specification	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	Floors should consists of wood, concrete or should be covered with ceramic tiles. If the floor is covered with synthetic material, the relative humidity must be at least 30 %.
Electrical fast Transient / Burst IEC 61000-4-4	± 2 kV for AC power line ± 1kV for input and output line 100 kHz	
Surges IEC 61000-4-5	± 1 kV differential mode ± 2kV common mode	
Voltage dips IEC 61000-4-11	0% U _T for ½ cycle (1 phase) 0% U _T for 1 cycle 70% U _T for 25/30 cycles ^a	
short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% U _T for 250/300 cycles ^a	
Power Frequency (50/60 Hz) Magnetic Field IEC 61000-4-8	30 A/m	
Note	U _T is the a.c. mains voltage prior to application of the test level.	
a	25/250 cycles at 50 Hz or 30/300 cycles at 60 Hz	

Table 2: Table 4, 5, 7, 8 EN 60601-1-2, Electromagnetic Immunity

MiniScreen PRO

Guidelines and Manufacturer's Declaration – Electromagnetic Immunity

The MiniScreen PRO is intended for an electromagnetic environment as specified below. The user of the MiniScreen PRO should make sure that is used in such an environment.

Immunity test	Test Specification	Electromagnetic Environment - Guidance
Conducted RF disturbance IEC 61000-4-6	3 V _{eff} 150kHz - 80MHz 6 V _{eff} ISM/Amateur bands 80 % AM / 1 kHz	Interferences may occur in the environment of devices, which are marked with the following symbol: 
Radiated RF disturbance IEC 61000-4-3	10 V/m 80MHz - 2,7 GHz 80 % AM / 1 kHz	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	380 - 390 MHz 27 V/m; PM 50%; 18 Hz	
	430 - 470 MHz 28 V/m; (FM ±5 kHz, 1 kHz sine) PM; 18 Hz	
	704 - 787 MHz 9 V/m; PM 50%; 217 Hz	
	800 - 960 MHz 28 V/m; PM 50%; 18 Hz	
	1700 - 1990 MHz 28 V/m; PM 50%; 217 Hz	
	2400 - 2570 MHz 28 V/m; PM 50%; 217 Hz	
	5100 - 5800 MHz 9 V/m; PM 50%; 217 Hz	

Table 3: Table 4, 5, 7, 8 EN 60601-1-2, Electromagnetic Immunity

MiniScreen PRO

Guidelines and Manufacturer's Declaration – Electromagnetic Immunity						
Test frequency MHz	Frequency band ^a MHz	Radio service ^a	Modulation ^b	Max. Performance W	Distance m	Immunity test level V/m
385	380 - 390	TETRA 400	Pulse Modulation ^b 18 Hz	1,8	0,3	27
450	430 - 470	GMRS 460, FRS 460	FM ^c ± 5 kHz Hub 1 kHz Sinus	2	0,3	28
710	704 - 787	LTE Band 13, 17	Pulse Modulation ^b 27 Hz	0,2	0,3	9
745						
780						
810	800 - 960	GSM 800/900 TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse Modulation ^b 18 Hz	2	0,3	28
870						
930						
1 720	1 700 - 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse Modulation ^b 217 Hz	2	0,3	28
1 845						
1 970						
2 450	2 400 - 2 570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse Modulation ^b 217 Hz	2	0,3	28
5 240	5 100 - 5 800	WLAN 802.11 a/n	Pulse Modulation ^b 217 Hz	0,2	0,3	9
5 500						
5 785						
NOTE: If necessary, to achieve immunity test levels, the distance between the transmit antenna and the MiniScreen PRO can be reduced to 1 m. The 1 m test distance is permitted according to IEC 61000-4-3.						
a	For some radio services, only the frequency for the radio link from the mobile communication device to the base station (en: uplink) has been included in the table. The carrier must be modulated with a square wave signal with 50% duty cycle.					
b	The carrier must be modulated with a square wave signal with 50% duty cycle.					
c	As an alternative to Frequency Modulation (FM), pulse modulation at 50% duty cycle at 18 Hz may be used, as this, if not the actual modulation, would be the worst case scenario					

Table 4: Table 9 EN 60601-1-2, Test specifications for the immunity of enclosures to high-frequency wireless communication devices

MiniScreen PRO

11. Index

Abdominal sensor.....	5
Accessories.....	2, 4
Assigned Purpose.....	2
Attaching the MiniScreen.....	7
Attaching the sensors.....	6
Battery charging.....	6
Charging the battery.....	7
Cleaning instructions.....	8
Disposal.....	4
EEG.....	6
EMG.....	6
EOG.....	6
Equipment.....	2
Finger sensor.....	5
Flow prong.....	5
Flow sensor.....	5
Functional testing.....	4
Information.....	3
Installing the software.....	9
Leg movement.....	6
Maintenance.....	9
Optical isolator.....	3
Ordering informations.....	11
PAP.....	2
PAP sensor.....	5
Position sensor.....	5
Preparing for measurement.....	6
Pressure pads.....	5
Safety-related information.....	3
Selecting a printer.....	9
Sensors.....	5
Sleep Xpert App.....	9
Starting the measurement.....	7
System requirements.....	2
Technical inspection.....	4
Technical specifications.....	12
Thermistor.....	5
Third-party accessories.....	4
Third-party software.....	4
Thorax sensor.....	5
Transport and storage conditions.....	9
Troubleshooting.....	9
Warranty.....	4
Whole system.....	2

Löwenstein Medical SE & Co. KG
Arzbacher Straße 80
56130 Bad Ems - Germany
Phone: +49 (0) 2603 9600-0
Fax: +49 (0) 2603 9600-40

CE 0197