

Health Products

AZEER
HEALTHCARE IS OUR MISSION



Healthcare Solutions



Best-Selling Products

- Specialized in providing medica-level products for every family in the world
- Annual sales more than 10 millions US dollars
- Top 3 in health industry of china



One M800 Series

Handheld Patient Monitor



Power saving



Adjustable alarms



300 groups of memory



Easy to carry

**Easier to carry
Easier to operate**

- Suitable for adult, pediatric and neonate
- Apply to hospital, clinic and home
- Monitoring parameter: SpO2,PR(Standard), Nellcor SpO2(Optional), ECG, RR(Optional)
- 2.4" color TFT display with auto rotation function
- Selectable layout modes including font and waveforms
- Dual work modes: continuous monitoring & spot check
- 3x AA batteries or rechargeable battery with charger stand(optional)
- Data can be transferred by communication software

M850 Handheld Patient Monitor



- Dual work modes
- Adjustable alarms
- 500 groups of memory
- Easy to carry

M860 Handheld Patient Monitor



**Easier to carry
Easier to operate**

- Dual work modes
- Adjustable alarms
- 500 groups of memory
- Easy to carry

- Suitable for adult, pediatric and neonate
- Apply to hospital, clinic and home
- Monitoring parameter: NIBP(Standard), SpO2, PR(Optional)
- 4.3"color TFT display with auto rotation function
- Specifically designed for a reliable and accurate blood pressure measurement
- Lithium battery support more than 20 hours continuous working
- Data can be transferred by communication software

**Easier to carry
Easier to operate**

- Suitable for adult, pediatric and neonate
- Apply to hospital, clinic and home
- Monitoring parameter: ECG, RR(Standard),SpO2, PR(Optional)
- 4.3"color TFT display with auto rotation function
- Selectable 3-lead and 5-lead for ECG measurement
- Selectable layout modes including font and waveforms
- Lithium battery support more than 28 hours continuous working
- Data can be transferred by communication software

M880 Handheld Patient Monitor



Dual work modes



Adjustable alarms



500 groups of memory



Easy to carry



Carry Bag
- M800



Case
- M800



Charge Stand
- M800



Lithium Battery
- M800



Trolley
- M850 M860 M880



Carry Bag
- M850 M860 M880



Case
- M850 M860 M880



Adult SpO2 Sensor
- M800 M850 M860 M880



Adult SpO2 Soft Reusable Sensor
- M800 M850 M860 M880



Pediatric SpO2 Reusable Sensor
- M800 M850 M860 M880



Neonatal SpO2 Soft Reusable Sensor
- M800 M850 M860 M880



Nellcor Adult SpO2 Sensor
- M800 M850 M860 M880

Easier to carry
Easier to operate

- Suitable for adult, pediatric and neonate
- Apply to hospital, clinic and home
- Monitoring parameter: CO2, RR(Standard), SpO2, PR(Optional)
- 4.3" color TFT display with auto rotation function
- Minimal warm-up time for EtCO2 and RR measurements
- CO2 water filter ensuring extended monitoring in a humid environment
- Selectable layout modes including font and waveforms
- More than 18 hours lithium battery
- Data can be transferred by communication software

Accessories - M800



Nellcor Neonatal SpO2 Sensor
- M800 M850 M860 M880



5-lead ECG Cable
- M850



3-lead ECG Cable
- M800 M850



ECG Adhesive Electrodes
- M800 M850



NIBP Adult Cuff
- M860



NIBP Large Adult Cuff
- M860



NIBP Pediatric Cuff
- M860



NIBP Neonatal Cuff
- M860



CO2 Tube
- M880



CO2 with 3-way Joint
- M880



CO2 Filter
- M880

Two M70 Series

Fingertip Pulse Oximeter

M70 fingertip Pulse Oximeter

M70A fingertip Pulse Oximeter

Care about your SpO2
Care about your health

Care about your SpO2
Care about your health



Medical level



Low power consumption



Medical level



Low power consumption



Front



Side

Key Features:

- Dual color OLED display
- Parameter: SpO2, PR, waveform and pulse bar
- 2 directions display
- Resistance to surrounding light and man-made light
- Available colors: Yellow, White
- Auto shut off without signal after 8s
- Suitable for pediatrics, adult and adolescent

Key Features:

- Cartoon design especially for children
- Dual color OLED display
- Parameter: SpO2, PR, waveform and pulse bar
- 2 directions display
- Resistance to surrounding light and man-made light
- Available colors: Orange, Pink
- Auto shut off without signal after 8s
- Suitable for pediatrics



Front



Side

M70B fingertip Pulse Oximeter

M70C fingertip Pulse Oximeter

Care about your SpO2
Care about your health



Medical level



Low power consumption



Front



Side

Key Features:

- OLED screen with 6 display modes
- 4 directions display
- Parameter: SpO2, PR, PI, waveform and pulse bar
- Resistance to surrounding light and man-made light
- Streamline design
- Auto shut off without signal after 8s
- Value signal includes Perfusion Index (PI)
- 5 adjustable brightness levels
- Suitable for pediatrics, adult and adolescent

Care about your SpO2
Care about your health



Medical level



Low power consumption



Easy to carry



Bluetooth (Optional)

Key Features:

- OLED screen with 6 display modes
- 4 directions display
- Colors optional: Blue, Yellow
- Parameter: SpO2, PR, waveform and pulse bar
- 5 adjustable brightness levels
- Resistance to surrounding light and man-made light
- Auto shut off without signal after 8s
- With Bluetooth and professional APP optional
- Suitable for pediatrics, adult and adolescent



Front



Side

M70D fingertip Pulse Oximeter



Care about your SpO2
Care about your health



Medical level



Low power consumption



Front



Side

Key Features:







- LCD display with backlight control
- Parameter: SpO2, PR and pulse bar
- Streamline design
- With soft silicon pad
- Resistance to surrounding light and man-made light
- Auto shut off without signal after 8s
- Suitable for pediatrics, adult and adolescent

Three WBP Series

Fingertip Pulse Oximeter

WBP 202 Blood Pressure Monitor



-  Bluetooth 4.0
-  One button operation
-  Alert function
-  Two users
-  Error detecting
-  Smart App







Simpler & Faster

- 2.8" LCD display
- All-in-one design
- Large data storage: 500 groups
- Intelligent pressurized
- Clinically accurate & fast reading
- Indication for wrong operation
- Various cuff options
- Operating time: > 24 hours at fully charge




WBP 302 Blood Pressure Monitor



-  Bluetooth 4.0
-  One button operation
-  Alert function
-  Two users
-  Error detecting
-  Smart App

Simpler & Faster

- Blood pressure, HR, Blood range display
- Large & eye protection backlight LCD screen
- Intelligent pressurized
- Clinically accurate & fast reading
- Indication for wrong operation
- Special design to place cuff
- 100 groups of memory
- Various cuff options



WF100 ultrasonic doppler fetal detector

Four WF 100

Ultrasonic doppler fetal heartbeat detector

Hearing your baby
from the beginning



Low radiation



Battery status indicator



Low power
consumption



Rechargeable and
equipped with adapter

- Compact and portable design Large color LCD display
- High sensitive probe
- Built-in speaker and interface of earphone
- Three working modes:
Real-time, Average and Manual counting
- Rechargeable Lithium battery

