



Circuits, ET T & LMA



Disposable Breathing Circuit
Product Code: BOT124000

Application	used together with anesthesia drug, respiratory machine, humidification device, and atomizer, to establish a respiratory channel for patients
Feature	1.Soft, flexible and air tight;
	2.All lengths and models available;
	3.Full sizes of connectors and extension tubes available;
	4.Suitable for breathing and anesthesia machine with different brand.



Model		Length (m) * length can be customized*							
Common Corrugated	Common Corrugated	1.0	1.2	1.4	1.5	1.6	1.8	2.0	
	W/ Water Trap	1.0	1.2	1.4	1.5	1.6	1.8	2.0	
	W/ Breathing Bag	1.0	1.2	1.4	1.5	1.6	1.8	2.0	
Smoothbore	Smoothbore	1.0	1.2	1.4	1.5	1.6	1.8	2.0	
	W/ Water Trap	1.0	1.2	1.4	1.5	1.6	1.8	2.0	
	W/ Breathing Bag	1.0	1.2	1.4	1.5	1.6	1.8	2.0	
Expandable	Expandable	1.0	1.2	1.4	1.5	1.6	1.8	2.0	
	W/ Water Trap	1.0	1.2	1.4	1.5	1.6	1.8	2.0	
	W/ Breathing Bag	1.0	1.2	1.4	1.5	1.6	1.8	2.0	
Co-axial	Co-axial	1.0	1.2	1.4	1.5	1.6	1.8	2.0	
	W/ Breathing Bag	1.0	1.2	1.4	1.5	1.6	1.8	2.0	



HME Filter

Product Code: BOT 111000

Application	used for filtering particles inside air ventilated by anesthesia machine and respiratory machine
Type	Straight and Angled
Model	Adult and Pediatric
Feature	1. Medical-grade high quality ABS material for outer shell;
	2. High-efficiency filtering and humidification medium;
	3. Good sealing, low breathing flow resistance and small dead space;
	4. Filtering efficiency more than 99.9%;
	5. Suitable for breathing and anesthesia machines with different brand;
	6. Easy to use.





Anesthesia Mask
Product Code: BOT 123000

Application	to establish and maintain an airway in general anesthesia
Feature	1. Made of clear and soft medical-grade material;
	2. Latex free;
	3. Air cushion ensures comfortable face fitting;
	4. Different size with different color code;
	5. 15mm-22mm standard connector.

Size	1#	2#	3#	4#	5#	6#
Connector	15mm	15mm	22mm	22mm	22mm	22mm

