



Oxygen therapy, Suction therapy, Air Oxygen Blenders



OXYGEN & SUCTION THERAPY

Flowmeters with Humidifiers | Medical Gas Regulators | Air/Oxygen Blenders | Suction Regulators | Collection Systems

Solutions for Life

SUCTION THERAPY

Suction Regulators

Solutions for Life

88XVR Series Suction Regulators

Mode Selector Switch

- Clearly marked indicator
- GENTEC® suction regulators are available in continuous only (882 Series) and continuous-intermittent (881 Series)
- 160 mmHg and 200 mmHg are available in pediatric, 300 mmHg is available in adult, and 760 mmHg is available in surgical
- Factory default settings for INT (on / off) cycles on the 881 Series is 8 seconds off, 16 seconds on. On/off settings can be adjusted

Control Knob

- Easy adjustment of vacuum pressure
- No Leakage
- Agion® Antimicrobial protection
 - Built-in long lasting antimicrobial protection
 - Based on naturally-occurring silver
 - Inhibits growth of microbes on surface of knobs and switches
 - Environmentally friendly technology
 - Broad Regulatory Approvals including EPA, FDA, NSF, USDA, Canada, European BPD, Swiss and EFSA

Fittings

- GENTEC® provides a wide range of fittings to meet customers' needs

Suction Trap

- Prevents contamination caused by overflow



Glow in the Dark Gauge

- 2-1/2" gauge, easy to read
- Color coded flow ranges
- Dual Scale mmHg and kPa
- Clockwise, counter-clockwise or digital gauge mode (features an analog scale in digital mode where a digital reading is on screen and an analog meter on the outer ring)



Clockwise Counter-clockwise Digital gauge

Face Cover

- High-quality ABS case body for long durability
- Impact-resistant internal components minimize damage if regulator is bumped or dropped
- Easy to clean, wipe down with sterilant

Interior Structure

- Few mechanical parts for easy maintenance
- Corrosion free
- Lightweight

881VR-300D-XX
Continuous / Intermittent
0-300 mmHg, digital gauge




















* Agion's silver antimicrobial has been certified for its material content, recyclability, and manufacturing characteristics. Cradle To Cradle Certified CM products meet established standards for human health, environmental health and recyclability. Cradle To Cradle Certified CM is a certification mark of MBDC.

SUCTION THERAPY

Solutions for Life

Suction Regulators

881VR Series Continuous / Intermittent				882VR Series Continuous-only			
Clockwise	 881VR-160 pediatric	 881VR-200 adult	 881VR-300 adult	 882VR-160 pediatric	 882VR-200 adult	 882VR-300 adult	 882VR-760 surgical
Counter-Clockwise	 881VR-160E pediatric		 881VR-300E adult	 882VR-160E adult		 882VR-300E adult	 882VR-760E surgical
Digital Gauge	 881VR-160D pediatric		 881VR-300D adult	 882VR-160D pediatric		 882VR-300D adult	 882VR-760D surgical
Mode Selection:							
OFF	√	√	√	√	√	√	√
CONT	√	√	√	√	√	√	√
INT	√	√	√	X	X	X	X
FULL	X	X	X	X	√	√	√
Regulated Vacuum Range	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~760 mmHg
Gauge Display Range	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~160 mmHg	0~200 mmHg	0~300 mmHg	0~760 mmHg
Gauge Accuracy	± 4%	± 4%	± 4%	± 4%	± 4%	± 4%	± 4%
Intermittent Time	Suction: 16 seconds, Stop: 8 seconds			/	/	/	/
Flow Rate	>50 l/min (CONT) 7-9 l/min (INT)	>55 l/min (CONT) 7-9 l/min (INT)	>60 l/min (CONT) 7-9 l/min (INT)	>50 l/min	>55 l/min	>60 l/min	>70 l/min
Weight	1 lb / 0.455 kg			0.81 lb / 0.37 kg			
Dimensions	6.3" H × 3.15" W × 3.94" D / 16 × 8 × 10 cm			6.3" H × 3.15" W × 3.94" D / 16 × 8 × 10 cm			

SUCTION THERAPY

Suction Regulators

Solutions for Life

88XVR Series Suction Regulators (Range 160)

 <p>Vacuum Connector</p> <p>Suction Connector</p>	 <p>No Adapter 1/8" NPT (F)</p>	 <p>OH: Omeda® 3111-VAC-NT2</p>	 <p>DH: DISS Hand-tight 3222-VAC-NT2</p>	 <p>DS: DISS Hex 3102-VAC-NT2</p>	 <p>CH: Chemetron® 3103-VAC-NT2</p>
 <p>No Adapter 1/8" NPT (F)</p>	<p>881VR-160-00-A 881VR-160E-00-A 881VR-160D-00-A 882VR-160-00-A 882VR-160E-00-A 882VR-160D-00-A 882VRM-160-00-A 882VRM-160E-00-A 882VRM-160D-00-A</p>	<p>881VR-160-OH-A 881VR-160E-OH-A 881VR-160D-OH-A 882VR-160-OH-A 882VR-160E-OH-A 882VR-160D-OH-A 882VRM-160-OH-A 882VRM-160E-OH-A 882VRM-160D-OH-A</p>	<p>881VR-160-DH-A 881VR-160E-DH-A 881VR-160D-DH-A 882VR-160-DH-A 882VR-160E-DH-A 882VR-160D-DH-A 882VRM-160-DH-A 882VRM-160E-DH-A 882VRM-160D-DH-A</p>	<p>881VR-160-DS-A 881VR-160E-DS-A 881VR-160D-DS-A 882VR-160-DS-A 882VR-160E-DS-A 882VR-160D-DS-A 882VRM-160-DS-A 882VRM-160E-DS-A 882VRM-160D-DS-A</p>	<p>881VR-160-CH-A 881VR-160E-CH-A 881VR-160D-CH-A 882VR-160-CH-A 882VR-160E-CH-A 882VR-160D-CH-A 882VRM-160-CH-A 882VRM-160E-CH-A 882VRM-160D-CH-A</p>
 <p>Suction Trap Assembly (complete)</p>	<p>881VR-160-00-T 881VR-160E-00-T 881VR-160D-00-T 882VR-160-00-T 882VR-160E-00-T 882VR-160D-00-T 882VRM-160-00-T 882VRM-160E-00-T 882VRM-160D-00-T</p>	<p>881VR-160-OH-T 881VR-160E-OH-T 881VR-160D-OH-T 882VR-160-OH-T 882VR-160E-OH-T 882VR-160D-OH-T 882VRM-160-OH-T 882VRM-160E-OH-T 882VRM-160D-OH-T</p>	<p>881VR-160-DH-T 881VR-160E-DH-T 881VR-160D-DH-T 882VR-160-DH-T 882VR-160E-DH-T 882VR-160D-DH-T 882VRM-160-DH-T 882VRM-160E-DH-T 882VRM-160D-DH-T</p>	<p>881VR-160-DS-T 881VR-160E-DS-T 881VR-160D-DS-T 882VR-160-DS-T 882VR-160E-DS-T 882VR-160D-DS-T 882VRM-160-DS-T 882VRM-160E-DS-T 882VRM-160D-DS-T</p>	<p>881VR-160-CH-T 881VR-160E-CH-T 881VR-160D-CH-T 882VR-160-CH-T 882VR-160E-CH-T 882VR-160D-CH-T 882VRM-160-CH-T 882VRM-160E-CH-T 882VRM-160D-CH-T</p>
 <p>Trap Fitting Only (No Trap)</p>	<p>881VR-160-00-B 881VR-160E-00-B 881VR-160D-00-B 882VR-160-00-B 882VR-160E-00-B 882VR-160D-00-B 882VRM-160-00-B 882VRM-160E-00-B 882VRM-160D-00-B</p>	<p>881VR-160-OH-B 881VR-160E-OH-B 881VR-160D-OH-B 882VR-160-OH-B 882VR-160E-OH-B 882VR-160D-OH-B 882VRM-160-OH-B 882VRM-160E-OH-B 882VRM-160D-OH-B</p>	<p>881VR-160-DH-B 881VR-160E-DH-B 881VR-160D-DH-B 882VR-160-DH-B 882VR-160E-DH-B 882VR-160D-DH-B 882VRM-160-DH-B 882VRM-160E-DH-B 882VRM-160D-DH-B</p>	<p>881VR-160-DS-B 881VR-160E-DS-B 881VR-160D-DS-B 882VR-160-DS-B 882VR-160E-DS-B 882VR-160D-DS-B 882VRM-160-DS-B 882VRM-160E-DS-B 882VRM-160D-DS-B</p>	<p>881VR-160-CH-B 881VR-160E-CH-B 881VR-160D-CH-B 882VR-160-CH-B 882VR-160E-CH-B 882VR-160D-CH-B 882VRM-160-CH-B 882VRM-160E-CH-B 882VRM-160D-CH-B</p>
 <p>DISS Male Vacuum</p>	<p>881VR-160-00-C 881VR-160E-00-C 881VR-160D-00-C 882VR-160-00-C 882VR-160E-00-C 882VR-160D-00-C 882VRM-160-00-C 882VRM-160E-00-C 882VRM-160D-00-C</p>	<p>881VR-160-OH-C 881VR-160E-OH-C 881VR-160D-OH-C 882VR-160-OH-C 882VR-160E-OH-C 882VR-160D-OH-C 882VRM-160-OH-C 882VRM-160E-OH-C 882VRM-160D-OH-C</p>	<p>881VR-160-DH-C 881VR-160E-DH-C 881VR-160D-DH-C 882VR-160-DH-C 882VR-160E-DH-C 882VR-160D-DH-C 882VRM-160-DH-C 882VRM-160E-DH-C 882VRM-160D-DH-C</p>	<p>881VR-160-DS-C 881VR-160E-DS-C 881VR-160D-DS-C 882VR-160-DS-C 882VR-160E-DS-C 882VR-160D-DS-C 882VRM-160-DS-C 882VRM-160E-DS-C 882VRM-160D-DS-C</p>	<p>881VR-160-CH-C 881VR-160E-CH-C 881VR-160D-CH-C 882VR-160-CH-C 882VR-160E-CH-C 882VR-160D-CH-C 882VRM-160-CH-C 882VRM-160E-CH-C 882VRM-160D-CH-C</p>
 <p>Tubing Nipple (5/16" Hose Barb)</p>	<p>881VR-160-00-D 881VR-160E-00-D 881VR-160D-00-D 882VR-160-00-D 882VR-160E-00-D 882VR-160D-00-D 882VRM-160-00-D 882VRM-160E-00-D 882VRM-160D-00-D</p>	<p>881VR-160-OH-D 881VR-160E-OH-D 881VR-160D-OH-D 882VR-160-OH-D 882VR-160E-OH-D 882VR-160D-OH-D 882VRM-160-OH-D 882VRM-160E-OH-D 882VRM-160D-OH-D</p>	<p>881VR-160-DH-D 881VR-160E-DH-D 881VR-160D-DH-D 882VR-160-DH-D 882VR-160E-DH-D 882VR-160D-DH-D 882VRM-160-DH-D 882VRM-160E-DH-D 882VRM-160D-DH-D</p>	<p>881VR-160-DS-D 881VR-160E-DS-D 881VR-160D-DS-D 882VR-160-DS-D 882VR-160E-DS-D 882VR-160D-DS-D 882VRM-160-DS-D 882VRM-160E-DS-D 882VRM-160D-DS-D</p>	<p>881VR-160-CH-D 881VR-160E-CH-D 881VR-160D-CH-D 882VR-160-CH-D 882VR-160E-CH-D 882VR-160D-CH-D 882VRM-160-CH-D 882VRM-160E-CH-D 882VRM-160D-CH-D</p>

					
PB: Puritan-Bennett® 3114-VAC-NT2	FS 3101-VAC-NT2	GS 3116-VAC-NT2	BS 3117-VAC-NT2	JIS 3108-VAC-NT2	SIS 3109-VAC-NT2
881VR-160-PB-A 881VR-160E-PB-A 881VR-160D-PB-A 882VR-160-PB-A 882VR-160E-PB-A 882VR-160D-PB-A 882VRM-160-PB-A 882VRM-160E-PB-A 882VRM-160D-PB-A	881VR-160-FS-A 881VR-160E-FS-A 881VR-160D-FS-A 882VR-160-FS-A 882VR-160E-FS-A 882VR-160D-FS-A 882VRM-160-FS-A 882VRM-160E-FS-A 882VRM-160D-FS-A	881VR-160-GS-A 881VR-160E-GS-A 881VR-160D-GS-A 882VR-160-GS-A 882VR-160E-GS-A 882VR-160D-GS-A 882VRM-160-GS-A 882VRM-160E-GS-A 882VRM-160D-GS-A	881VR-160-BS-A 881VR-160E-BS-A 881VR-160D-BS-A 882VR-160-BS-A 882VR-160E-BS-A 882VR-160D-BS-A 882VRM-160-BS-A 882VRM-160E-BS-A 882VRM-160D-BS-A	881VR-160-JIS-A 881VR-160E-JIS-A 881VR-160D-JIS-A 882VR-160-JIS-A 882VR-160E-JIS-A 882VR-160D-JIS-A 882VRM-160-JIS-A 882VRM-160E-JIS-A 882VRM-160D-JIS-A	881VR-160-SIS-A 881VR-160E-SIS-A 881VR-160D-SIS-A 882VR-160-SIS-A 882VR-160E-SIS-A 882VR-160D-SIS-A 882VRM-160-SIS-A 882VRM-160E-SIS-A 882VRM-160D-SIS-A
881VR-160-PB-T 881VR-160E-PB-T 881VR-160D-PB-T 882VR-160-PB-T 882VR-160E-PB-T 882VR-160D-PB-T 882VRM-160-PB-T 882VRM-160E-PB-T 882VRM-160D-PB-T	881VR-160-FS-T 881VR-160E-FS-T 881VR-160D-FS-T 882VR-160-FS-T 882VR-160E-FS-T 882VR-160D-FS-T 882VRM-160-FS-T 882VRM-160E-FS-T 882VRM-160D-FS-T	881VR-160-GS-T 881VR-160E-GS-T 881VR-160D-GS-T 882VR-160-GS-T 882VR-160E-GS-T 882VR-160D-GS-T 882VRM-160-GS-T 882VRM-160E-GS-T 882VRM-160D-GS-T	881VR-160-BS-T 881VR-160E-BS-T 881VR-160D-BS-T 882VR-160-BS-T 882VR-160E-BS-T 882VR-160D-BS-T 882VRM-160-BS-T 882VRM-160E-BS-T 882VRM-160D-BS-T	881VR-160-JIS-T 881VR-160E-JIS-T 881VR-160D-JIS-T 882VR-160-JIS-T 882VR-160E-JIS-T 882VR-160D-JIS-T 882VRM-160-JIS-T 882VRM-160E-JIS-T 882VRM-160D-JIS-T	881VR-160-SIS-T 881VR-160E-SIS-T 881VR-160D-SIS-T 882VR-160-SIS-T 882VR-160E-SIS-T 882VR-160D-SIS-T 882VRM-160-SIS-T 882VRM-160E-SIS-T 882VRM-160D-SIS-T
881VR-160-PB-B 881VR-160E-PB-B 881VR-160D-PB-B 882VR-160-PB-B 882VR-160E-PB-B 882VR-160D-PB-B 882VRM-160-PB-B 882VRM-160E-PB-B 882VRM-160D-PB-B	881VR-160-FS-B 881VR-160E-FS-B 881VR-160D-FS-B 882VR-160-FS-B 882VR-160E-FS-B 882VR-160D-FS-B 882VRM-160-FS-B 882VRM-160E-FS-B 882VRM-160D-FS-B	881VR-160-GS-B 881VR-160E-GS-B 881VR-160D-GS-B 882VR-160-GS-B 882VR-160E-GS-B 882VR-160D-GS-B 882VRM-160-GS-B 882VRM-160E-GS-B 882VRM-160D-GS-B	881VR-160-BS-B 881VR-160E-BS-B 881VR-160D-BS-B 882VR-160-BS-B 882VR-160E-BS-B 882VR-160D-BS-B 882VRM-160-BS-B 882VRM-160E-BS-B 882VRM-160D-BS-B	881VR-160-JIS-B 881VR-160E-JIS-B 881VR-160D-JIS-B 882VR-160-JIS-B 882VR-160E-JIS-B 882VR-160D-JIS-B 882VRM-160-JIS-B 882VRM-160E-JIS-B 882VRM-160D-JIS-B	881VR-160-SIS-B 881VR-160E-SIS-B 881VR-160D-SIS-B 882VR-160-SIS-B 882VR-160E-SIS-B 882VR-160D-SIS-B 882VRM-160-SIS-B 882VRM-160E-SIS-B 882VRM-160D-SIS-B
881VR-160-PB-C 881VR-160E-PB-C 881VR-160D-PB-C 882VR-160-PB-C 882VR-160E-PB-C 882VR-160D-PB-C 882VRM-160-PB-C 882VRM-160E-PB-C 882VRM-160D-PB-C	881VR-160-FS-C 881VR-160E-FS-C 881VR-160D-FS-C 882VR-160-FS-C 882VR-160E-FS-C 882VR-160D-FS-C 882VRM-160-FS-C 882VRM-160E-FS-C 882VRM-160D-FS-C	881VR-160-GS-C 881VR-160E-GS-C 881VR-160D-GS-C 882VR-160-GS-C 882VR-160E-GS-C 882VR-160D-GS-C 882VRM-160-GS-C 882VRM-160E-GS-C 882VRM-160D-GS-C	881VR-160-BS-C 881VR-160E-BS-C 881VR-160D-BS-C 882VR-160-BS-C 882VR-160E-BS-C 882VR-160D-BS-C 882VRM-160-BS-C 882VRM-160E-BS-C 882VRM-160D-BS-C	881VR-160-JIS-C 881VR-160E-JIS-C 881VR-160D-JIS-C 882VR-160-JIS-C 882VR-160E-JIS-C 882VR-160D-JIS-C 882VRM-160-JIS-C 882VRM-160E-JIS-C 882VRM-160D-JIS-C	881VR-160-SIS-C 881VR-160E-SIS-C 881VR-160D-SIS-C 882VR-160-SIS-C 882VR-160E-SIS-C 882VR-160D-SIS-C 882VRM-160-SIS-C 882VRM-160E-SIS-C 882VRM-160D-SIS-C
881VR-160-PB-D 881VR-160E-PB-D 881VR-160D-PB-D 882VR-160-PB-D 882VR-160E-PB-D 882VR-160D-PB-D 882VRM-160-PB-D 882VRM-160E-PB-D 882VRM-160D-PB-D	881VR-160-FS-D 881VR-160E-FS-D 881VR-160D-FS-D 882VR-160-FS-D 882VR-160E-FS-D 882VR-160D-FS-D 882VRM-160-FS-D 882VRM-160E-FS-D 882VRM-160D-FS-D	881VR-160-GS-D 881VR-160E-GS-D 881VR-160D-GS-D 882VR-160-GS-D 882VR-160E-GS-D 882VR-160D-GS-D 882VRM-160-GS-D 882VRM-160E-GS-D 882VRM-160D-GS-D	881VR-160-BS-D 881VR-160E-BS-D 881VR-160D-BS-D 882VR-160-BS-D 882VR-160E-BS-D 882VR-160D-BS-D 882VRM-160-BS-D 882VRM-160E-BS-D 882VRM-160D-BS-D	881VR-160-JIS-D 881VR-160E-JIS-D 881VR-160D-JIS-D 882VR-160-JIS-D 882VR-160E-JIS-D 882VR-160D-JIS-D 882VRM-160-JIS-D 882VRM-160E-JIS-D 882VRM-160D-JIS-D	881VR-160-SIS-D 881VR-160E-SIS-D 881VR-160D-SIS-D 882VR-160-SIS-D 882VR-160E-SIS-D 882VR-160D-SIS-D 882VRM-160-SIS-D 882VRM-160E-SIS-D 882VRM-160D-SIS-D

see
previous
page for
complete
table

SUCTION THERAPY

Suction Regulators

Solutions for Life

88XVR Series Suction Regulators (Range 200)

 <p>Vacuum Connector</p> <p>Suction Connector</p>	 <p>No Adapter 1/8" NPT (F)</p>	 <p>OH: Omeda® 3111-VAC-NT2</p>	 <p>DH: DISS Hand-tight 3222-VAC-NT2</p>	 <p>DS: DISS Hex 3102-VAC-NT2</p>	 <p>CH: Chemetron® 3103-VAC-NT2</p>
 <p>No Adapter 1/8" NPT (F)</p>	<p>881VR-200-00-A 882VR-200-00-A 882VRM-200-00-A</p>	<p>881VR-200-OH-A 882VR-200-OH-A 882VRM-200-OH-A</p>	<p>881VR-200-DH-A 882VR-200-DH-A 882VRM-200-DH-A</p>	<p>881VR-200-DS-A 882VR-200-DS-A 882VRM-200-DS-A</p>	<p>881VR-200-CH-A 882VR-200-CH-A 882VRM-200-CH-A</p>
 <p>Suction Trap Assembly (complete)</p>	<p>881VR-200-00-T 882VR-200-00-T 882VRM-200-00-T</p>	<p>881VR-200-OH-T 882VR-200-OH-T 882VRM-200-OH-T</p>	<p>881VR-200-DH-T 882VR-200-DH-T 882VRM-200-DH-T</p>	<p>881VR-200-DS-T 882VR-200-DS-T 882VRM-200-DS-T</p>	<p>881VR-200-CH-T 882VR-200-CH-T 882VRM-200-CH-T</p>
 <p>Trap Fitting Only (No Trap)</p>	<p>881VR-200-00-B 882VR-200-00-B 882VRM-200-00-B</p>	<p>881VR-200-OH-B 882VR-200-OH-B 882VRM-200-OH-B</p>	<p>881VR-200-DH-B 882VR-200-DH-B 882VRM-200-DH-B</p>	<p>881VR-200-DS-B 882VR-200-DS-B 882VRM-200-DS-B</p>	<p>881VR-200-CH-B 882VR-200-CH-B 882VRM-200-CH-B</p>
 <p>DISS Male Vacuum</p>	<p>881VR-200-00-C 882VR-200-00-C 882VRM-200-00-C</p>	<p>881VR-200-OH-C 882VR-200-OH-C 882VRM-200-OH-C</p>	<p>881VR-200-DH-C 882VR-200-DH-C 882VRM-200-DH-C</p>	<p>881VR-200-DS-C 882VR-200-DS-C 882VRM-200-DS-C</p>	<p>881VR-200-CH-C 882VR-200-CH-C 882VRM-200-CH-C</p>
 <p>Tubing Nipple (5/16" Hose Barb)</p>	<p>881VR-200-00-D 882VR-200-00-D 882VRM-200-00-D</p>	<p>881VR-200-OH-D 882VR-200-OH-D 882VRM-200-OH-D</p>	<p>881VR-200-DH-D 882VR-200-DH-D 882VRM-200-DH-D</p>	<p>881VR-200-DS-D 882VR-200-DS-D 882VRM-200-DS-D</p>	<p>881VR-200-CH-D 882VR-200-CH-D 882VRM-200-CH-D</p>

 <p>PB: Puritan-Bennett® 3114-VAC-NT2</p>	 <p>FS 3101-VAC-NT2</p>	 <p>GS 3116-VAC-NT2</p>	 <p>BS 3117-VAC-NT2</p>	 <p>JIS 3108-VAC-NT2</p>	 <p>SIS 3109-VAC-NT2</p>
<p>881VR-200-PB-A 882VR-200-PB-A 882VRM-200-PB-A</p>	<p>881VR-200-FS-A 882VR-200-FS-A 882VRM-200-FS-A</p>	<p>881VR-200-GS-A 882VR-200-GS-A 882VRM-200-GS-A</p>	<p>881VR-200-BS-A 882VR-200-BS-A 882VRM-200-BS-A</p>	<p>881VR-200-JIS-A 882VR-200-JIS-A 882VRM-200-JIS-A</p>	<p>881VR-200-SIS-A 882VR-200-SIS-A 882VRM-200-SIS-A</p>
<p>881VR-200-PB-T 882VR-200-PB-T 882VRM-200-PB-T</p>	<p>881VR-200-FS-T 882VR-200-FS-T 882VRM-200-FS-T</p>	<p>881VR-200-GS-T 882VR-200-GS-T 882VRM-200-GS-T</p>	<p>881VR-200-BS-T 882VR-200-BS-T 882VRM-200-BS-T</p>	<p>881VR-200-JIS-T 882VR-200-JIS-T 882VRM-200-JIS-T</p>	<p>881VR-200-SIS-T 882VR-200-SIS-T 882VRM-200-SIS-T</p>
<p>881VR-200-PB-B 882VR-200-PB-B 882VRM-200-PB-B</p>	<p>881VR-200-FS-B 882VR-200-FS-B 882VRM-200-FS-B</p>	<p>881VR-200-GS-B 882VR-200-GS-B 882VRM-200-GS-B</p>	<p>881VR-200-BS-B 882VR-200-BS-B 882VRM-200-BS-B</p>	<p>881VR-200-JIS-B 882VR-200-JIS-B 882VRM-200-JIS-B</p>	<p>881VR-200-SIS-B 882VR-200-SIS-B 882VRM-200-SIS-B</p>
<p>881VR-200-PB-C 882VR-200-PB-C 882VRM-200-PB-C</p>	<p>881VR-200-FS-C 882VR-200-FS-C 882VRM-200-FS-C</p>	<p>881VR-200-GS-C 882VR-200-GS-C 882VRM-200-GS-C</p>	<p>881VR-200-BS-C 882VR-200-BS-C 882VRM-200-BS-C</p>	<p>881VR-200-JIS-C 882VR-200-JIS-C 882VRM-200-JIS-C</p>	<p>881VR-200-SIS-C 882VR-200-SIS-C 882VRM-200-SIS-C</p>
<p>881VR-200-PB-D 882VR-200-PB-D 882VRM-200-PB-D</p>	<p>881VR-200-FS-D 882VR-200-FS-D 882VRM-200-FS-D</p>	<p>881VR-200-GS-D 882VR-200-GS-D 882VRM-200-GS-D</p>	<p>881VR-200-BS-D 882VR-200-BS-D 882VRM-200-BS-D</p>	<p>881VR-200-JIS-D 882VR-200-JIS-D 882VRM-200-JIS-D</p>	<p>881VR-200-SIS-D 882VR-200-SIS-D 882VRM-200-SIS-D</p>

see
previous
page for
complete
table


SUCTION THERAPY

Suction Regulators

Solutions for Life

88XVR Series Suction Regulators (Range 300)

 <p>Vacuum Connector</p> <p>Suction Connector</p>	 <p>No Adapter 1/8" NPT (F)</p>	 <p>OH: Omeda® 3111-VAC-NT2</p>	 <p>DH: DISS Hand-tight 3222-VAC-NT2</p>	 <p>DS: DISS Hex 3102-VAC-NT2</p>	 <p>CH: Chemetron® 3103-VAC-NT2</p>
 <p>No Adapter 1/8" NPT (F)</p>	<p>881VR-300-00-A 881VR-300E-00-A 881VR-300D-00-A 882VR-300-00-A 882VR-300E-00-A 882VR-300D-00-A 882VRM-300-00-A 882VRM-300E-00-A 882VRM-300D-00-A</p>	<p>881VR-300-OH-A 881VR-300E-OH-A 881VR-300D-OH-A 882VR-300-OH-A 882VR-300E-OH-A 882VR-300D-OH-A 882VRM-300-OH-A 882VRM-300E-OH-A 882VRM-300D-OH-A</p>	<p>881VR-300-DH-A 881VR-300E-DH-A 881VR-300D-DH-A 882VR-300-DH-A 882VR-300E-DH-A 882VR-300D-DH-A 882VRM-300-DH-A 882VRM-300E-DH-A 882VRM-300D-DH-A</p>	<p>881VR-300-DS-A 881VR-300E-DS-A 881VR-300D-DS-A 882VR-300-DS-A 882VR-300E-DS-A 882VR-300D-DS-A 882VRM-300-DS-A 882VRM-300E-DS-A 882VRM-300D-DS-A</p>	<p>881VR-300-CH-A 881VR-300E-CH-A 881VR-300D-CH-A 882VR-300-CH-A 882VR-300E-CH-A 882VR-300D-CH-A 882VRM-300-CH-A 882VRM-300E-CH-A 882VRM-300D-CH-A</p>
 <p>Suction Trap Assembly (complete)</p>	<p>881VR-300-00-T 881VR-300E-00-T 881VR-300D-00-T 882VR-300-00-T 882VR-300E-00-T 882VR-300D-00-T 882VRM-300-00-T 882VRM-300E-00-T 882VRM-300D-00-T</p>	<p>881VR-300-OH-T 881VR-300E-OH-T 881VR-300D-OH-T 882VR-300-OH-T 882VR-300E-OH-T 882VR-300D-OH-T 882VRM-300-OH-T 882VRM-300E-OH-T 882VRM-300D-OH-T</p>	<p>881VR-300-DH-T 881VR-300E-DH-T 881VR-300D-DH-T 882VR-300-DH-T 882VR-300E-DH-T 882VR-300D-DH-T 882VRM-300-DH-T 882VRM-300E-DH-T 882VRM-300D-DH-T</p>	<p>881VR-300-DS-T 881VR-300E-DS-T 881VR-300D-DS-T 882VR-300-DS-T 882VR-300E-DS-T 882VR-300D-DS-T 882VRM-300-DS-T 882VRM-300E-DS-T 882VRM-300D-DS-T</p>	<p>881VR-300-CH-T 881VR-300E-CH-T 881VR-300D-CH-T 882VR-300-CH-T 882VR-300E-CH-T 882VR-300D-CH-T 882VRM-300-CH-T 882VRM-300E-CH-T 882VRM-300D-CH-T</p>
 <p>Trap Fitting Only (No Trap)</p>	<p>881VR-300-00-B 881VR-300E-00-B 881VR-300D-00-B 882VR-300-00-B 882VR-300E-00-B 882VR-300D-00-B 882VRM-300-00-B 882VRM-300E-00-B 882VRM-300D-00-B</p>	<p>881VR-300-OH-B 881VR-300E-OH-B 881VR-300D-OH-B 882VR-300-OH-B 882VR-300E-OH-B 882VR-300D-OH-B 882VRM-300-OH-B 882VRM-300E-OH-B 882VRM-300D-OH-B</p>	<p>881VR-300-DH-B 881VR-300E-DH-B 881VR-300D-DH-B 882VR-300-DH-B 882VR-300E-DH-B 882VR-300D-DH-B 882VRM-300-DH-B 882VRM-300E-DH-B 882VRM-300D-DH-B</p>	<p>881VR-300-DS-B 881VR-300E-DS-B 881VR-300D-DS-B 882VR-300-DS-B 882VR-300E-DS-B 882VR-300D-DS-B 882VRM-300-DS-B 882VRM-300E-DS-B 882VRM-300D-DS-B</p>	<p>881VR-300-CH-B 881VR-300E-CH-B 881VR-300D-CH-B 882VR-300-CH-B 882VR-300E-CH-B 882VR-300D-CH-B 882VRM-300-CH-B 882VRM-300E-CH-B 882VRM-300D-CH-B</p>
 <p>DISS Male Vacuum</p>	<p>881VR-300-00-C 881VR-300E-00-C 881VR-300D-00-C 882VR-300-00-C 882VR-300E-00-C 882VR-300D-00-C 882VRM-300-00-C 882VRM-300E-00-C 882VRM-300D-00-C</p>	<p>881VR-300-OH-C 881VR-300E-OH-C 881VR-300D-OH-C 882VR-300-OH-C 882VR-300E-OH-C 882VR-300D-OH-C 882VRM-300-OH-C 882VRM-300E-OH-C 882VRM-300D-OH-C</p>	<p>881VR-300-DH-C 881VR-300E-DH-C 881VR-300D-DH-C 882VR-300-DH-C 882VR-300E-DH-C 882VR-300D-DH-C 882VRM-300-DH-C 882VRM-300E-DH-C 882VRM-300D-DH-C</p>	<p>881VR-300-DS-C 881VR-300E-DS-C 881VR-300D-DS-C 882VR-300-DS-C 882VR-300E-DS-C 882VR-300D-DS-C 882VRM-300-DS-C 882VRM-300E-DS-C 882VRM-300D-DS-C</p>	<p>881VR-300-CH-C 881VR-300E-CH-C 881VR-300D-CH-C 882VR-300-CH-C 882VR-300E-CH-C 882VR-300D-CH-C 882VRM-300-CH-C 882VRM-300E-CH-C 882VRM-300D-CH-C</p>
 <p>Tubing Nipple (5/16" Hose Barb)</p>	<p>881VR-300-00-D 881VR-300E-00-D 881VR-300D-00-D 882VR-300-00-D 882VR-300E-00-D 882VR-300D-00-D 882VRM-300-00-D 882VRM-300E-00-D 882VRM-300D-00-D</p>	<p>881VR-300-OH-D 881VR-300E-OH-D 881VR-300D-OH-D 882VR-300-OH-D 882VR-300E-OH-D 882VR-300D-OH-D 882VRM-300-OH-D 882VRM-300E-OH-D 882VRM-300D-OH-D</p>	<p>881VR-300-DH-D 881VR-300E-DH-D 881VR-300D-DH-D 882VR-300-DH-D 882VR-300E-DH-D 882VR-300D-DH-D 882VRM-300-DH-D 882VRM-300E-DH-D 882VRM-300D-DH-D</p>	<p>881VR-300-DS-D 881VR-300E-DS-D 881VR-300D-DS-D 882VR-300-DS-D 882VR-300E-DS-D 882VR-300D-DS-D 882VRM-300-DS-D 882VRM-300E-DS-D 882VRM-300D-DS-D</p>	<p>881VR-300-CH-D 881VR-300E-CH-D 881VR-300D-CH-D 882VR-300-CH-D 882VR-300E-CH-D 882VR-300D-CH-D 882VRM-300-CH-D 882VRM-300E-CH-D 882VRM-300D-CH-D</p>

					
PB: Puritan-Bennett® 3114-VAC-NT2	FS 3101-VAC-NT2	GS 3116-VAC-NT2	BS 3117-VAC-NT2	JIS 3108-VAC-NT2	SIS 3109-VAC-NT2
881VR-300-PB-A 881VR-300E-PB-A 881VR-300D-PB-A 882VR-300-PB-A 882VR-300E-PB-A 882VR-300D-PB-A 882VRM-300-PB-A 882VRM-300E-PB-A 882VRM-300D-PB-A	881VR-300-FS-A 881VR-300E-FS-A 881VR-300D-FS-A 882VR-300-FS-A 882VR-300E-FS-A 882VR-300D-FS-A 882VRM-300-FS-A 882VRM-300E-FS-A 882VRM-300D-FS-A	881VR-300-GS-A 881VR-300E-GS-A 881VR-300D-GS-A 882VR-300-GS-A 882VR-300E-GS-A 882VR-300D-GS-A 882VRM-300-GS-A 882VRM-300E-GS-A 882VRM-300D-GS-A	881VR-300-BS-A 881VR-300E-BS-A 881VR-300D-BS-A 882VR-300-BS-A 882VR-300E-BS-A 882VR-300D-BS-A 882VRM-300-BS-A 882VRM-300E-BS-A 882VRM-300D-BS-A	881VR-300-JIS-A 881VR-300E-JIS-A 881VR-300D-JIS-A 882VR-300-JIS-A 882VR-300E-JIS-A 882VR-300D-JIS-A 882VRM-300-JIS-A 882VRM-300E-JIS-A 882VRM-300D-JIS-A	881VR-300-SIS-A 881VR-300E-SIS-A 881VR-300D-SIS-A 882VR-300-SIS-A 882VR-300E-SIS-A 882VR-300D-SIS-A 882VRM-300-SIS-A 882VRM-300E-SIS-A 882VRM-300D-SIS-A
881VR-300-PB-T 881VR-300E-PB-T 881VR-300D-PB-T 882VR-300-PB-T 882VR-300E-PB-T 882VR-300D-PB-T 882VRM-300-PB-T 882VRM-300E-PB-T 882VRM-300D-PB-T	881VR-300-FS-T 881VR-300E-FS-T 881VR-300D-FS-T 882VR-300-FS-T 882VR-300E-FS-T 882VR-300D-FS-T 882VRM-300-FS-T 882VRM-300E-FS-T 882VRM-300D-FS-T	881VR-300-GS-T 881VR-300E-GS-T 881VR-300D-GS-T 882VR-300-GS-T 882VR-300E-GS-T 882VR-300D-GS-T 882VRM-300-GS-T 882VRM-300E-GS-T 882VRM-300D-GS-T	881VR-300-BS-T 881VR-300E-BS-T 881VR-300D-BS-T 882VR-300-BS-T 882VR-300E-BS-T 882VR-300D-BS-T 882VRM-300-BS-T 882VRM-300E-BS-T 882VRM-300D-BS-T	881VR-300-JIS-T 881VR-300E-JIS-T 881VR-300D-JIS-T 882VR-300-JIS-T 882VR-300E-JIS-T 882VR-300D-JIS-T 882VRM-300-JIS-T 882VRM-300E-JIS-T 882VRM-300D-JIS-T	881VR-300-SIS-T 881VR-300E-SIS-T 881VR-300D-SIS-T 882VR-300-SIS-T 882VR-300E-SIS-T 882VR-300D-SIS-T 882VRM-300-SIS-T 882VRM-300E-SIS-T 882VRM-300D-SIS-T
881VR-300-PB-B 881VR-300E-PB-B 881VR-300D-PB-B 882VR-300-PB-B 882VR-300E-PB-B 882VR-300D-PB-B 882VRM-300-PB-B 882VRM-300E-PB-B 882VRM-300D-PB-B	881VR-300-FS-B 881VR-300E-FS-B 881VR-300D-FS-B 882VR-300-FS-B 882VR-300E-FS-B 882VR-300D-FS-B 882VRM-300-FS-B 882VRM-300E-FS-B 882VRM-300D-FS-B	881VR-300-GS-B 881VR-300E-GS-B 881VR-300D-GS-B 882VR-300-GS-B 882VR-300E-GS-B 882VR-300D-GS-B 882VRM-300-GS-B 882VRM-300E-GS-B 882VRM-300D-GS-B	881VR-300-BS-B 881VR-300E-BS-B 881VR-300D-BS-B 882VR-300-BS-B 882VR-300E-BS-B 882VR-300D-BS-B 882VRM-300-BS-B 882VRM-300E-BS-B 882VRM-300D-BS-B	881VR-300-JIS-B 881VR-300E-JIS-B 881VR-300D-JIS-B 882VR-300-JIS-B 882VR-300E-JIS-B 882VR-300D-JIS-B 882VRM-300-JIS-B 882VRM-300E-JIS-B 882VRM-300D-JIS-B	881VR-300-SIS-B 881VR-300E-SIS-B 881VR-300D-SIS-B 882VR-300-SIS-B 882VR-300E-SIS-B 882VR-300D-SIS-B 882VRM-300-SIS-B 882VRM-300E-SIS-B 882VRM-300D-SIS-B
881VR-300-PB-C 881VR-300E-PB-C 881VR-300D-PB-C 882VR-300-PB-C 882VR-300E-PB-C 882VR-300D-PB-C 882VRM-300-PB-C 882VRM-300E-PB-C 882VRM-300D-PB-C	881VR-300-FS-C 881VR-300E-FS-C 881VR-300D-FS-C 882VR-300-FS-C 882VR-300E-FS-C 882VR-300D-FS-C 882VRM-300-FS-C 882VRM-300E-FS-C 882VRM-300D-FS-C	881VR-300-GS-C 881VR-300E-GS-C 881VR-300D-GS-C 882VR-300-GS-C 882VR-300E-GS-C 882VR-300D-GS-C 882VRM-300-GS-C 882VRM-300E-GS-C 882VRM-300D-GS-C	881VR-300-BS-C 881VR-300E-BS-C 881VR-300D-BS-C 882VR-300-BS-C 882VR-300E-BS-C 882VR-300D-BS-C 882VRM-300-BS-C 882VRM-300E-BS-C 882VRM-300D-BS-C	881VR-300-JIS-C 881VR-300E-JIS-C 881VR-300D-JIS-C 882VR-300-JIS-C 882VR-300E-JIS-C 882VR-300D-JIS-C 882VRM-300-JIS-C 882VRM-300E-JIS-C 882VRM-300D-JIS-C	881VR-300-SIS-C 881VR-300E-SIS-C 881VR-300D-SIS-C 882VR-300-SIS-C 882VR-300E-SIS-C 882VR-300D-SIS-C 882VRM-300-SIS-C 882VRM-300E-SIS-C 882VRM-300D-SIS-C
881VR-300-PB-D 881VR-300E-PB-D 881VR-300D-PB-D 882VR-300-PB-D 882VR-300E-PB-D 882VR-300D-PB-D 882VRM-300-PB-D 882VRM-300E-PB-D 882VRM-300D-PB-D	881VR-300-FS-D 881VR-300E-FS-D 881VR-300D-FS-D 882VR-300-FS-D 882VR-300E-FS-D 882VR-300D-FS-D 882VRM-300-FS-D 882VRM-300E-FS-D 882VRM-300D-FS-D	881VR-300-GS-D 881VR-300E-GS-D 881VR-300D-GS-D 882VR-300-GS-D 882VR-300E-GS-D 882VR-300D-GS-D 882VRM-300-GS-D 882VRM-300E-GS-D 882VRM-300D-GS-D	881VR-300-BS-D 881VR-300E-BS-D 881VR-300D-BS-D 882VR-300-BS-D 882VR-300E-BS-D 882VR-300D-BS-D 882VRM-300-BS-D 882VRM-300E-BS-D 882VRM-300D-BS-D	881VR-300-JIS-D 881VR-300E-JIS-D 881VR-300D-JIS-D 882VR-300-JIS-D 882VR-300E-JIS-D 882VR-300D-JIS-D 882VRM-300-JIS-D 882VRM-300E-JIS-D 882VRM-300D-JIS-D	881VR-300-SIS-D 881VR-300E-SIS-D 881VR-300D-SIS-D 882VR-300-SIS-D 882VR-300E-SIS-D 882VR-300D-SIS-D 882VRM-300-SIS-D 882VRM-300E-SIS-D 882VRM-300D-SIS-D

see
previous
page for
complete
table

SUCTION THERAPY

Suction Regulators

Solutions for Life

88XVR Series Suction Regulators (Range 760)

 <p>Vacuum Connector</p> <p>Suction Connector</p>	 <p>No Adapter 1/8" NPT (F)</p>	 <p>OH: Omeda® 3111-VAC-NT2</p>	 <p>DH: DISS Hand-tight 3222-VAC-NT2</p>	 <p>DS: DISS Hex 3102-VAC-NT2</p>	 <p>CH: Chemetron® 3103-VAC-NT2</p>
 <p>No Adapter 1/8" NPT (F)</p>	<p>882VR-760-00-A 882VR-760E-00-A 882VR-760D-00-A 882VRM-760-00-A 882VRM-760E-00-A 882VRM-760D-00-A</p>	<p>882VR-760-OH-A 882VR-760E-OH-A 882VR-760D-OH-A 882VRM-760-OH-A 882VRM-760E-OH-A 882VRM-760D-OH-A</p>	<p>882VR-760-DH-A 882VR-760E-DH-A 882VR-760D-DH-A 882VRM-760-DH-A 882VRM-760E-DH-A 882VRM-760D-DH-A</p>	<p>882VR-760-DS-A 882VR-760E-DS-A 882VR-760D-DS-A 882VRM-760-DS-A 882VRM-760E-DS-A 882VRM-760D-DS-A</p>	<p>882VR-760-CH-A 882VR-760E-CH-A 882VR-760D-CH-A 882VRM-760-CH-A 882VRM-760E-CH-A 882VRM-760D-CH-A</p>
 <p>Suction Trap Assembly (complete)</p>	<p>882VR-760-00-T 882VR-760E-00-T 882VR-760D-00-T 882VRM-760-00-T 882VRM-760E-00-T 882VRM-760D-00-T</p>	<p>882VR-760-OH-T 882VR-760E-OH-T 882VR-760D-OH-T 882VRM-760-OH-T 882VRM-760E-OH-T 882VRM-760D-OH-T</p>	<p>882VR-760-DH-T 882VR-760E-DH-T 882VR-760D-DH-T 882VRM-760-DH-T 882VRM-760E-DH-T 882VRM-760D-DH-T</p>	<p>882VR-760-DS-T 882VR-760E-DS-T 882VR-760D-DS-T 882VRM-760-DS-T 882VRM-760E-DS-T 882VRM-760D-DS-T</p>	<p>882VR-760-CH-T 882VR-760E-CH-T 882VR-760D-CH-T 882VRM-760-CH-T 882VRM-760E-CH-T 882VRM-760D-CH-T</p>
 <p>Trap Fitting Only (No Trap)</p>	<p>882VR-760-00-B 882VR-760E-00-B 882VR-760D-00-B 882VRM-760-00-B 882VRM-760E-00-B 882VRM-760D-00-B</p>	<p>882VR-760-OH-B 882VR-760E-OH-B 882VR-760D-OH-B 882VRM-760-OH-B 882VRM-760E-OH-B 882VRM-760D-OH-B</p>	<p>882VR-760-DH-B 882VR-760E-DH-B 882VR-760D-DH-B 882VRM-760-DH-B 882VRM-760E-DH-B 882VRM-760D-DH-B</p>	<p>882VR-760-DS-B 882VR-760E-DS-B 882VR-760D-DS-B 882VRM-760-DS-B 882VRM-760E-DS-B 882VRM-760D-DS-B</p>	<p>882VR-760-CH-B 882VR-760E-CH-B 882VR-760D-CH-B 882VRM-760-CH-B 882VRM-760E-CH-B 882VRM-760D-CH-B</p>
 <p>DISS Male Vacuum</p>	<p>882VR-760-00-C 882VR-760E-00-C 882VR-760D-00-C 882VRM-760-00-C 882VRM-760E-00-C 882VRM-760D-00-C</p>	<p>882VR-760-OH-C 882VR-760E-OH-C 882VR-760D-OH-C 882VRM-760-OH-C 882VRM-760E-OH-C 882VRM-760D-OH-C</p>	<p>882VR-760-DH-C 882VR-760E-DH-C 882VR-760D-DH-C 882VRM-760-DH-C 882VRM-760E-DH-C 882VRM-760D-DH-C</p>	<p>882VR-760-DS-C 882VR-760E-DS-C 882VR-760D-DS-C 882VRM-760-DS-C 882VRM-760E-DS-C 882VRM-760D-DS-C</p>	<p>882VR-760-CH-C 882VR-760E-CH-C 882VR-760D-CH-C 882VRM-760-CH-C 882VRM-760E-CH-C 882VRM-760D-CH-C</p>
 <p>Tubing Nipple (5/16" Hose Barb)</p>	<p>882VR-760-00-D 882VR-760E-00-D 882VR-760D-00-D 882VRM-760-00-D 882VRM-760E-00-D 882VRM-760D-00-D</p>	<p>882VR-760-OH-D 882VR-760E-OH-D 882VR-760D-OH-D 882VRM-760-OH-D 882VRM-760E-OH-D 882VRM-760D-OH-D</p>	<p>882VR-760-DH-D 882VR-760E-DH-D 882VR-760D-DH-D 882VRM-760-DH-D 882VRM-760E-DH-D 882VRM-760D-DH-D</p>	<p>882VR-760-DS-D 882VR-760E-DS-D 882VR-760D-DS-D 882VRM-760-DS-D 882VRM-760E-DS-D 882VRM-760D-DS-D</p>	<p>882VR-760-CH-D 882VR-760E-CH-D 882VR-760D-CH-D 882VRM-760-CH-D 882VRM-760E-CH-D 882VRM-760D-CH-D</p>

 <p>PB: Puritan-Bennett® 3114-VAC-NT2</p>	 <p>FS 3101-VAC-NT2</p>	 <p>GS 3116-VAC-NT2</p>	 <p>BS 3117-VAC-NT2</p>	 <p>JIS 3108-VAC-NT2</p>	 <p>SIS 3109-VAC-NT2</p>
882VR-760-PB-A 882VR-760E-PB-A 882VR-760D-PB-A 882VRM-760-PB-A 882VRM-760E-PB-A 882VRM-760D-PB-A	882VR-760-FS-A 882VR-760E-FS-A 882VR-760D-FS-A 882VRM-760-FS-A 882VRM-760E-FS-A 882VRM-760D-FS-A	882VR-760-GS-A 882VR-760E-GS-A 882VR-760D-GS-A 882VRM-760-GS-A 882VRM-760E-GS-A 882VRM-760D-GS-A	882VR-760-BS-A 882VR-760E-BS-A 882VR-760D-BS-A 882VRM-760-BS-A 882VRM-760E-BS-A 882VRM-760D-BS-A	882VR-760-JIS-A 882VR-760E-JIS-A 882VR-760D-JIS-A 882VRM-760-JIS-A 882VRM-760E-JIS-A 882VRM-760D-JIS-A	882VR-760-SIS-A 882VR-760E-SIS-A 882VR-760D-SIS-A 882VRM-760-SIS-A 882VRM-760E-SIS-A 882VRM-760D-SIS-A
882VR-760-PB-T 882VR-760E-PB-T 882VR-760D-PB-T 882VRM-760-PB-T 882VRM-760E-PB-T 882VRM-760D-PB-T	882VR-760-FS-T 882VR-760E-FS-T 882VR-760D-FS-T 882VRM-760-FS-T 882VRM-760E-FS-T 882VRM-760D-FS-T	882VR-760-GS-T 882VR-760E-GS-T 882VR-760D-GS-T 882VRM-760-GS-T 882VRM-760E-GS-T 882VRM-760D-GS-T	882VR-760-BS-T 882VR-760E-BS-T 882VR-760D-BS-T 882VRM-760-BS-T 882VRM-760E-BS-T 882VRM-760D-BS-T	882VR-760-JIS-T 882VR-760E-JIS-T 882VR-760D-JIS-T 882VRM-760-JIS-T 882VRM-760E-JIS-T 882VRM-760D-JIS-T	882VR-760-SIS-T 882VR-760E-SIS-T 882VR-760D-SIS-T 882VRM-760-SIS-T 882VRM-760E-SIS-T 882VRM-760D-SIS-T
882VR-760-PB-B 882VR-760E-PB-B 882VR-760D-PB-B 882VRM-760-PB-B 882VRM-760E-PB-B 882VRM-760D-PB-B	882VR-760-FS-B 882VR-760E-FS-B 882VR-760D-FS-B 882VRM-760-FS-B 882VRM-760E-FS-B 882VRM-760D-FS-B	882VR-760-GS-B 882VR-760E-GS-B 882VR-760D-GS-B 882VRM-760-GS-B 882VRM-760E-GS-B 882VRM-760D-GS-B	882VR-760-BS-B 882VR-760E-BS-B 882VR-760D-BS-B 882VRM-760-BS-B 882VRM-760E-BS-B 882VRM-760D-BS-B	882VR-760-JIS-B 882VR-760E-JIS-B 882VR-760D-JIS-B 882VRM-760-JIS-B 882VRM-760E-JIS-B 882VRM-760D-JIS-B	882VR-760-SIS-B 882VR-760E-SIS-B 882VR-760D-SIS-B 882VRM-760-SIS-B 882VRM-760E-SIS-B 882VRM-760D-SIS-B
882VR-760-PB-C 882VR-760E-PB-C 882VR-760D-PB-C 882VRM-760-PB-C 882VRM-760E-PB-C 882VRM-760D-PB-C	882VR-760-FS-C 882VR-760E-FS-C 882VR-760D-FS-C 882VRM-760-FS-C 882VRM-760E-FS-C 882VRM-760D-FS-C	882VR-760-GS-C 882VR-760E-GS-C 882VR-760D-GS-C 882VRM-760-GS-C 882VRM-760E-GS-C 882VRM-760D-GS-C	882VR-760-BS-C 882VR-760E-BS-C 882VR-760D-BS-C 882VRM-760-BS-C 882VRM-760E-BS-C 882VRM-760D-BS-C	882VR-760-JIS-C 882VR-760E-JIS-C 882VR-760D-JIS-C 882VRM-760-JIS-C 882VRM-760E-JIS-C 882VRM-760D-JIS-C	882VR-760-SIS-C 882VR-760E-SIS-C 882VR-760D-SIS-C 882VRM-760-SIS-C 882VRM-760E-SIS-C 882VRM-760D-SIS-C
882VR-760-PB-D 882VR-760E-PB-D 882VR-760D-PB-D 882VRM-760-PB-D 882VRM-760E-PB-D 882VRM-760D-PB-D	882VR-760-FS-D 882VR-760E-FS-D 882VR-760D-FS-D 882VRM-760-FS-D 882VRM-760E-FS-D 882VRM-760D-FS-D	882VR-760-GS-D 882VR-760E-GS-D 882VR-760D-GS-D 882VRM-760-GS-D 882VRM-760E-GS-D 882VRM-760D-GS-D	882VR-760-BS-D 882VR-760E-BS-D 882VR-760D-BS-D 882VRM-760-BS-D 882VRM-760E-BS-D 882VRM-760D-BS-D	882VR-760-JIS-D 882VR-760E-JIS-D 882VR-760D-JIS-D 882VRM-760-JIS-D 882VRM-760E-JIS-D 882VRM-760D-JIS-D	882VR-760-SIS-D 882VR-760E-SIS-D 882VR-760D-SIS-D 882VRM-760-SIS-D 882VRM-760E-SIS-D 882VRM-760D-SIS-D

see
previous
page for
complete
table

SUCTION THERAPY

Suction Regulator

Solutions for Life

880VS Series Venturi Suction Regulator



Glow in the Dark Gauge

- 2" gauge, easy to read
- Color coded flow ranges
- Dual Scale mmHg and kPa

Safety Device

- Positive pressure relief valve provides safe use
- ISO 10079-3 Compliant

Control Knob

- Easy adjustment of vacuum pressure
- No Leakage

Suction Trap

- Prevents contamination caused by overflow

Type of Adapter

- None: Without adapter
- OH: Ohmeda® adapter
- DS: DISS hex adapter
- DH: DISS Handwheel adapter
- CH: Chemetron® adapter
- PB: Puritan-Bennett® adapter
- FS: French (NF S 90-116) adapter
- GS: German (DIN 13260-2) adapter
- BS: British (BS5682-1998) adapter
- JIS: Japanese Style adapter
- SIS: Australian (AS2896) adapter

Venturi Suction

- Delivers steady and high vacuum flow from a medical air or oxygen gas source
- Control Range: -1 kPa ~ -80 kPa

Suction Connector

- GENTEC® provides a wide range of fittings to meet customers' needs



880VS-760-00-A-T

Technical Specification

Dimension (H x W x D)	90 x 90 x 60 mm
Weight	490 g
Range	-1 kPa ~ -80 kPa
Maximum Suction	> 25 L / min
Source of Pressure	300 ~ 500 kPa
Vacuum Control	Needle valve control venturi
Vacuum Gauge Range	0 ~ -100 kPa / 0 ~ 760 mmHg

Ordering Information

Ex: 880VS	- 760	- 00	- A	- T
Series	Measuring Range	Vacuum Connection	Driving Gas	Suction Connection
880VS	760: 0~760 mmHg	00: No adapter, 1/8" NPT (F) OH: 3111U (Ohmeda) adapter DH: 3222U (DISS hand-tight) adapter DS: 3102 (DISS Hex) adapter CH: 3103 (Chemetron) adapter PB: 3104 (Puritan-Bennett) adapter FS: 3105 (FS) adapter GS: 3116 (GS) adapter BS: 3117 (BS) adapter JIS: 3108 (JIS) adapter SIS: 3109 (SIS) adapter	A: Medical Air O: Oxygen	00: No adapter, 1/8" NPT (F) T: With suction trap assembly (complete) TD: With suction trap (complete DISS Female nut) TH: With suction trap (complete DISS female vacuum hand-tight) H: With tubing nipple (5/16" hose barb) D: With DISS male vacuum N: With NIST adapter F: With trap fitting only (no trap)

SUCTION THERAPY

Suction Regulator

Solutions for Life

880VS Series Venturi Suction Regulator

 Vacuum Connector  Suction Connector	 No Adapter No Adapter, 1/8" NPT (F)	 OH: Omeda® 3101-OXY-NT2 3111-AIR-NT2	 DH: DISS Hand-tight 3222-OXY-NT2 3222-AIR-NT2	 DS: DISS Hex 3102-OXY-NT2 3102-AIR-NT2	 CH: Chemetron® 3103-OXY-NT2 3103-AIR-NT2
 No Adapter No Adapter, 1/8" NPT (F)	880VS-760-00-A-00 880VS-760-00-O-00	880VS-760-OH-A-00 880VS-760-OH-O-00	880VS-760-DH-A-00 880VS-760-DH-O-00	880VS-760-DS-A-00 880VS-760-DS-O-00	880VS-760-CH-A-00 880VS-760-CH-O-00
 Suction Trap Assembly(complete)	880VS-760-00-A-T 880VS-760-00-O-T	880VS-760-OH-A-T 880VS-760-OH-O-T	880VS-760-DH-A-T 880VS-760-DH-O-T	880VS-760-DS-A-T 880VS-760-DS-O-T	880VS-760-CH-A-T 880VS-760-CH-O-T
 Suction Trap (with DISS female nut & nipple vacuum fitting)	880VS-760-00-A-TD 880VS-760-00-O-TD	880VS-760-OH-A-TD 880VS-760-OH-O-TD	880VS-760-DH-A-TD 880VS-760-DH-O-TD	880VS-760-DS-A-TD 880VS-760-DS-O-TD	880VS-760-CH-A-TD 880VS-760-CH-O-TD
 Suction Trap (with DISS female vacuum hand-tight)	880VS-760-00-A-TH 880VS-760-00-O-TH	880VS-760-OH-A-TH 880VS-760-OH-O-TH	880VS-760-DH-A-TH 880VS-760-DH-O-TH	880VS-760-DS-A-TH 880VS-760-DS-O-TH	880VS-760-CH-A-TH 880VS-760-CH-O-TH
 Tubing Nipple (5/16" Hose Barb)	880VS-760-00-A-H 880VS-760-00-O-H	880VS-760-OH-A-H 880VS-760-OH-O-H	880VS-760-DH-A-H 880VS-760-DH-O-H	880VS-760-DS-A-H 880VS-760-DS-O-H	880VS-760-CH-A-H 880VS-760-CH-O-H
 DISS Male Vacuum	880VS-760-00-A-D 880VS-760-00-O-D	880VS-760-OH-A-D 880VS-760-OH-O-D	880VS-760-DH-A-D 880VS-760-DH-O-D	880VS-760-DS-A-D 880VS-760-DS-O-D	880VS-760-CH-A-D 880VS-760-CH-O-D
 NIST vacuum fitting	880VS-760-00-A-N 880VS-760-00-O-N	880VS-760-OH-A-N 880VS-760-OH-O-N	880VS-760-DH-A-N 880VS-760-DH-O-N	880VS-760-DS-A-N 880VS-760-DS-O-N	880VS-760-CH-A-N 880VS-760-CH-O-N
 Trap Fitting Only (No Trap)	880VS-760-00-A-F 880VS-760-00-O-F	880VS-760-OH-A-F 880VS-760-OH-O-F	880VS-760-DH-A-F 880VS-760-DH-O-F	880VS-760-DS-A-F 880VS-760-DS-O-F	880VS-760-CH-A-F 880VS-760-CH-O-F

 PB: Puritan-Bennett® 3104-OXY-NT2 3114-AIR-NT2	 FS 3101-OXY-NT2 3101-AIR-NT2	 GS 3116-OXY-NT2 3116-AIR-NT2	 BS 3117-OXY-NT2 3117-AIR-NT2	 JIS 3108-OXY-NT2 3108-AIR-NT2	 SIS 3109-OXY-NT2 3109-AIR-NT2
880VS-760-PB-A-00 880VS-760-PB-O-00	880VS-760-FS-A-00 880VS-760-FS-O-00	880VS-760-GS-A-00 880VS-760-GS-O-00	880VS-760-BS-A-00 880VS-760-BS-O-00	880VS-760-JIS-A-00 880VS-760-JIS-O-00	880VS-760-SIS-A-00 880VS-760-SIS-O-00
880VS-760-PB-A-T 880VS-760-PB-O-T	880VS-760-FS-A-T 880VS-760-FS-O-T	880VS-760-GS-A-T 880VS-760-GS-O-T	880VS-760-BS-A-T 880VS-760-BS-O-T	880VS-760-JIS-A-T 880VS-760-JIS-O-T	880VS-760-SIS-A-T 880VS-760-SIS-O-T
880VS-760-PB-A-TD 880VS-760-PB-O-TD	880VS-760-FS-A-TD 880VS-760-FS-O-TD	880VS-760-GS-A-TD 880VS-760-GS-O-TD	880VS-760-BS-A-TD 880VS-760-BS-O-TD	880VS-760-JIS-A-TD 880VS-760-JIS-O-TD	880VS-760-SIS-A-TD 880VS-760-SIS-O-TD
880VS-760-PB-A-TH 880VS-760-PB-O-TH	880VS-760-FS-A-TH 880VS-760-FS-O-TH	880VS-760-GS-A-TH 880VS-760-GS-O-TH	880VS-760-BS-A-TH 880VS-760-BS-O-TH	880VS-760-JIS-A-TH 880VS-760-JIS-O-TH	880VS-760-SIS-A-TH 880VS-760-SIS-O-TH
880VS-760-PB-A-H 880VS-760-PB-O-H	880VS-760-FS-A-H 880VS-760-FS-O-H	880VS-760-GS-A-H 880VS-760-GS-O-H	880VS-760-BS-A-H 880VS-760-BS-O-H	880VS-760-JIS-A-H 880VS-760-JIS-O-H	880VS-760-SIS-A-H 880VS-760-SIS-O-H
880VS-760-PB-A-D 880VS-760-PB-O-D	880VS-760-FS-A-D 880VS-760-FS-O-D	880VS-760-GS-A-D 880VS-760-GS-O-D	880VS-760-BS-A-D 880VS-760-BS-O-D	880VS-760-JIS-A-D 880VS-760-JIS-O-D	880VS-760-SIS-A-D 880VS-760-SIS-O-D
880VS-760-PB-A-N 880VS-760-PB-O-N	880VS-760-FS-A-N 880VS-760-FS-O-N	880VS-760-GS-A-N 880VS-760-GS-O-N	880VS-760-BS-A-N 880VS-760-BS-O-N	880VS-760-JIS-A-N 880VS-760-JIS-O-N	880VS-760-SIS-A-N 880VS-760-SIS-O-N
880VS-760-PB-A-F 880VS-760-PB-O-F	880VS-760-FS-A-F 880VS-760-FS-O-F	880VS-760-GS-A-F 880VS-760-GS-O-F	880VS-760-BS-A-F 880VS-760-BS-O-F	880VS-760-JIS-A-F 880VS-760-JIS-O-F	880VS-760-SIS-A-F 880VS-760-SIS-O-F

see
previous
page for
complete
table

SUCTION THERAPY

Accessories

Solutions for Life

Suction Traps

Suction traps prevent contamination caused by overflow

880VT-0	880VT-1	880VT-2	880VT-3
 <p data-bbox="261 963 428 1008">with no fitting on top (1/8" NPT female)</p>	 <p data-bbox="532 963 651 1008">with male trap fitting only</p>	 <p data-bbox="732 963 915 1008">with DISS female nut & nipple vacuum fitting</p>	 <p data-bbox="976 963 1127 1008">with DISS female vacuum hand-tight</p>

880VT-4	880VT-5	880VT-6	880VT-7	880VT-8
 <p data-bbox="272 1780 423 1824">with tubing nipple (5/16" hose barb)</p>	 <p data-bbox="483 1780 699 1824">with complete trap fitting assembly (male & female)</p>	 <p data-bbox="743 1759 911 1824">with complete DISS Female nut & nipple vacuum fitting</p>	 <p data-bbox="971 1759 1138 1824">with complete DISS female vacuum hand-tight</p>	 <p data-bbox="1208 1780 1338 1824">with small filter, disposable</p>

Suction Gauges

<p>Clockwise</p>	 <p>GR25PM-V160 (0 ~ 160 mmHg)</p>	 <p>GR25PM-V200 (0 ~ 200 mmHg)</p>	 <p>GR25PM-V300 (0 ~ 300 mmHg)</p>	 <p>GR25PM-V760 (0 ~ 760 mmHg)</p>
<p>Counter-Clockwise</p>	 <p>GR25PM-V160-E (0 ~ 160 mmHg)</p>	 <p>GR25PM-V300-E (0 ~ 300 mmHg)</p>	 <p>GR25PM-V760-E (0 ~ 760 mmHg)</p>	

Venturi Vacuum Modules

Features

- Designed to work with vacuum regulator
- Small and compact
- With positive pressure relief valve
- With silencer
- Maximum vacuum: 0 ~ 600 mmHg
- Maximum air free flow: 30 L / min
- Noise: < 60 dBa



880VE-X



SUCTION THERAPY

Medical Ejector Suction System

Solutions for Life

Venturi Vacuum System



Ordering Information

Ex: 881VR - 300	E	- VG	A	- 00	- T	- 1L	* 2	
Series	Measuring Range	Pointer Rotate Mode	Suction Generator	Driving Gas	Vacuum Connection	Suction Connection	Canister Volumn	Canister Quantity
881VR	160: 0~160 mmHg	None:	VG: With Venturi	A: Air	00: No adapter	T: Suction trap assembly (complete)	1L: 1L canister	1: one canister
882VR	200: 0~200 mmHg 300: 0~300 mmHg 760: 0~760 mmHg (only for 882VR)	Clockwise, (USA) E: Counter-clockwise (ISO)	Vacuum module None: Without Venturi Vacuum module	O: Oxygen	1/8" NPT (F) FDS: DISS hose barb (USA) FNS: NIST hose barb (ISO)	None: Without suction trap	1.5L: 1.5L canister 2.5L: 2.5L canister	2: two canisters

Standard suction hose length is 0.5 m, for other length, please contact GENTEC.

Disposable Bacteria Filters

- Prevents contamination and fluid from entering the vacuum regulator and vacuum systems
- PVDF special filter material, can effectively filter the virus and bacteria
- Effective filtration area: 32 cm²
- Maximum operating temperature: 140 °F (60 °C)
- Maximum operating pressure: -14.5 ~ 14.5 psi
- Filters down to 0.22 micron




F11-HB



F11-NT2/HB



Dimensions



F11-HB

F11-NT2/HB

Ordering Information

Model	Description
F11-NT2/HB-0.45	Hydrophobic, Inlet: 1/4" Hose Barb, Outlet: 1/8" NPT, Filter Accuracy: 0.45µm
F11-NT2/HB-0.22	Hydrophobic, Inlet: 1/4" Hose Barb, Outlet: 1/8" NPT, Filter Accuracy: 0.22µm
F11-HB-0.45	Hydrophobic, Inlet and Outlet: 1/4" Hose Barb, Filter Accuracy: 0.45µm
F11-HB-0.22	Hydrophobic, Inlet and Outlet: 1/4" Hose Barb, Filter Accuracy: 0.22µm



Genstar Technologies Company Inc.

4525 Edison Avenue, Chino, CA 91710, USA

Tel: 909-606-2726

Fax: 909-606-6485

www.gentechealthcare.com

